

# An Introduction to HPLC Instrumentation

## Part 1 Agilent HPLC Webinar Series

Lori Sandford, Applications Scientist  
Agilent Technologies, Inc.  
Wood Dale, IL



# Topics for Today's Seminar

Agilent 1260/1290 Specifications  
Instrument Flowpath

The Pump

- General Hygiene of Solvents
- Pump Design
- Starting up the Instrument/Prime/Purge
- Pumping Parameters
- Compositional Accuracy

The Autosampler

- Injection Mechanism

The Column Compartment Settings and  
Proper Use

HPLC Detectors



- UV Detection
  - UV Settings
  - Sensitivity

Definitions of Basic Parameters for  
Optimization

- Delay Volume
- Dispersion

# Typical HPLC System


Detector  
Column  
Compartment  
Autosampler  
Pump &  
Vacuum  
Degasser



# Typical HPLC System

## Follow the Flowpath

- Where are the moving parts?
- Where can blockages occur?
- Where are the consumable parts?
- Where can leaks occur ?




***What can I do to eliminate, reduce or anticipate potential problems with the LC ?***

# Filters and Bottle necks for blockages


## Reduce Downtime by Understanding Common Problems

- Solvent inlet filters in solvent bottles
  - glass: 20um – replace if needed!
  - SST: 12-14um – replace, opt sonicate
- Heat exchanger (bent? Connection?)
- PTFE frit in the purge valve at outlet of the quaternary pump head: replace
- Binary pump only: inline filter at outlet check valve
- Troubleshoot: Disconnect other modules behind pump



# The HPLC Stack

## Common Consumable Parts



Typical setup with **0.17mm ID green capillaries**

### Solvent bottles :

01018-60025 Solvent inlet filter **12-14um SS**

3150-0944 Glass filter, solvent inlet, 40um

5043-1218 A-Line Stay Safe cap (GL45), 2-port

9301-1450 Solvent bottle (GL45) amber, 1 L

9301-1421 Solvent bottle (GL45) clear, 1 L with cap

### Pump (1260 Quat with integrated degasser):

G1311B Agilent 1260 Infinity Quaternary Pump  
PTFE Frit (pack of 5) 01018-22707

### Autosampler - Column Compartment:

G1367E Agilent 1260 Infinity High Performance Autosampler

G1329B Agilent 1260 Standard Autosampler

### Column Compartment - Column:

G1330B Agilent 1260 Infinity Thermostatted Column Compartment

5067-5965 A-Line quick-connect LC fitting

### Detector - waste:

G4212B 1260 Infinity Diode-Array Detector

G4212B #030 1260 Max-Light Hi-sensitivity cell

G314F Variable Wavelength Detector

# Pumps

- General Hygiene and Housekeeping
- Preparing Pumps
- High and Low Pressure Mixing Pumps
- Compositional Accuracy
- Compressibility of Solvents
- Seal Wash




# Performance Characteristics of the Pump

## Important Characteristics

### Common to isocratic and gradient pumps

- Flow accuracy
- Flow precision
- Pressure pulsation




## Influence on...

- Retention time and peak area precision (system to system)
- Retention time and peak area precision (within one system)
- Baseline noise

### Common to gradient pumps only

- Delay volume in low and high pressure mixing
- Compositional accuracy
- Composition precision



- Gradient shape and precision
- Retention time and peak area precision (system to system)
- Retention time and peak area precision (within one system)



# Reduce LC problems by eliminating most common sources of flow blockage

- Filter or use HPLC grade solvents
- Particles leading to blockage can come from sources located both outside and inside the LC system:
  - Solvent, buffer
  - Microbial growth in solvent reservoirs
  - The Sample
  - Wear of LC components – piston seals, autosampler valve, etc.
- Debris will either be captured on the column frit or in-line filter
- Use 0.45  $\mu\text{m}$  filter for mobile phase components
- Replace aqueous solvents every two days
- Avoid exposure to direct sunlight
- Preventative Maintenance is the Key!
- Keep a log of maintenance to prevent downtime



# Examples: Used / Unused Filters


Glass filters: 3150 - 0944



Stainless Steel Filters: 01018 – 60025  
(less volume, no Na+ ions)




# Clean or Replace the Solvent Inlet Filters



Clogged inlet filter led to cavitation which in one line causing pH change which changed selectivity

# Agilent 1260 Infinity Quaternary Pump


## Low-pressure mixing (LPM) principle



Mixing by low-pressure proportioning valve before the pump head

# Agilent 1260 Infinity Binary Pump

## High-pressure mixing (HPM) principle



Combination and mixing of mobile phases after the pump head

# When to use purge, prime, condition?

## **Purge**

Change solvents

When pump is refilled with new/different mobile phase the purge valves allows both pump heads (binary pump) to be connected to waste at the same time

## **Prime**

When the pump is dry

When Purge and Condition still show exhausted pressure ripple

## **Condition**

When first starting up for the day or after changing solvents


When pump pressure ripple or composition ripple is too high (mixing noise) air bubble is hidden in pump head (listen)

best once a day to condition for smooth operation




# Solvent Compressibility

- Solvent compressibility must be accounted for, especially when running at UHPLC pressures >9000 psi
- Failure to account for compressibility will result in different flow rates at higher pressures



# Compositional Accuracy



Select the correct solvent type in your pump parameters



*Shallow gradient at low organic concentration. In this case the trace in red cannot perform the steps at low composition. This system should not be used for shallow gradients (<1% delta B / minute)*



# Many Pumps Have Seal Wash Option and an In-Line Filter



Pump head


Peristaltic Pump

- Protects your pistons and seals from excessive wear
- Recommended when using aqueous buffer or salt solutions  $> 0.1M$
- 10% isopropanol recommended as pump seal wash solvent
- Change solvent weekly (date on bottle)

# Pump and Degasser Maintenance

- Clean the degasser lines by flushing with isopropanol.
- When using buffers, flush with water, then with isopropanol.
- Check for air bubbles in outlet lines.
- Be aware of the possibility of microbial growth in aqueous phases
  - flush aqueous lines weekly with IPA
  - Don't top off water/solvents, use fresh bottles that have been cleaned with solvent and dried
- Check for solvent compatibility and flush with appropriate solvents
- Unused channels should be left in isopropanol.

# Pump Head – Main Components



# Autosampler – the brain behind workflow automation logistics

- Sampler Settings
  - Fixed Loop vs Flow-Through Autosamplers
  - Principle of Operation



# Performance characteristics of the Autosampler

## Characteristic

- Injection volume precision
- Wide linearity
- Minimum carry over
- Wide dynamic injection volume




## Influences

- Precision of peak area/height
- Accuracy of peak area/height (when using different injection volumes)
- Precision of peak area/height
- Versatility, application range

# Fixed Loop Autosampler

Sample is drawn or pushed into a sample loop and this loop is switched into the sample flow-path

- To achieve highest reproducibility you must usually overload the sample loop. Partial filling of the loop is possible but reproducibility of this technique is poor.
- To inject smaller or larger volumes with a fixed-loop autosampler, you must typically change the sample loop.




# Flow-Through Autosampler


## Most Agilent Autosamplers

The mobile phase passes through the sample loop as well through the metering device and injection needle.


- The major advantage of a flow-through design is more flexibility in terms of the injection volume range
- Flowing through the needle during the run helps to eliminate carryover




# Principle of Operation



**Prior to Injection**  
**Valve in Mainpass Position**



**Draw Sample**  
**Valve in Bypass Position**



**Injection and Run**  
**Valve in Mainpass Position**



# Thermostatted Column Compartment

## Important performance characteristics

Excellent temperature accuracy



Excellent temperature precision



## Influence on...

- Elution order
- Peak identification
  
- Elution order
- Retention time precision
- Peak identification


# Column Compartment

- Temperature Settings
- Precolumn Heating
  - Needed for analyses done at non-ambient condition
  - Columns with smaller IDs (<4.6mm) and shorter length (<250mm) will see greater temperature disparities from front to back of column due to insufficient mobile phase heating


High Efficiency Heat Exchangers



# Effect of Temperature on Separation



# HPLC Detector Usage



# HPLC Detector Characteristics

Detector Type	Sensitivity	Selectivity	Useful % of Compounds	Advantages
VWD	ng/pg	-	80	Low cost
DAD	ng/pg	++	80	Peak confirmation
Fluorescence	pg/fg	++	10	High sensitivity
Electro-chemical	pg/fg	+	>20	High sensitivity
Conductivity	ng/pg	-	10	Ion chromatography
Refractive Index	µg/ng	-	100	Universal response
Mass Spectrometer	ng/pg	++	<100	<b>MW structural information</b>

# UV-Based Detectors

Fixed or Variable Detector


- One Wavelength or programmable Wavelength

Diode Array Detector

- Multiwavelength
- Spectra (Library, Peak Purity)



# Variable Wavelength (VWD) UV/VIS Detector Optical System




**Figure 1** Optical Path of the Variable Wavelength Detector



# Diode Array Detector (DAD) UV/VIS Optical System


**Figure 59**      **Optical System of the Detector**







# Detector – Baseline noise

## Type and geometry of flow cell

4.6 mm id column + 10 mm path length 13  $\mu$ L volume = 

Short 2.1 mm id column + 10 mm path length 13  $\mu$ L volume = 

Short 2.1 mm id column + Long path length (10 mm, 0.5  $\mu$ L) Low light transmission. => high noise = 

# Data Rate Setting

Setup Method

Binary Pump Quat. Pump2 HiP Sampler HiP Sampler Injector Program Column Comp. DAD DAD2 Fraction Collector Instrument Curves

DAD (G4212A)

Signals

	Use Signal	Wave length	Band width	Reference Wavelength	Reference Bandwidth	
Signal A	<input checked="" type="checkbox"/>	254.0	4.0	<input checked="" type="checkbox"/> 360.0	100.0	nm
Signal B	<input checked="" type="checkbox"/>	210.0	4.0	<input checked="" type="checkbox"/> 360.0	100.0	nm
Signal C	<input checked="" type="checkbox"/>	214.0	4.0	<input checked="" type="checkbox"/> 360.0	100.0	nm
Signal D	<input checked="" type="checkbox"/>	230.0	4.0	<input checked="" type="checkbox"/> 360.0	100.0	nm
Signal E	<input checked="" type="checkbox"/>	260.0	4.0	<input checked="" type="checkbox"/> 360.0	100.0	nm
Signal F	<input type="checkbox"/>	273.0	4.0	<input checked="" type="checkbox"/> 360.0	100.0	nm
Signal G	<input type="checkbox"/>	280.0	4.0	<input checked="" type="checkbox"/> 360.0	100.0	nm
Signal H	<input type="checkbox"/>	250.0	4.0	<input checked="" type="checkbox"/> 360.0	100.0	nm

Peakwidth

>0.10 min (2.0 s response time) (2.5 Hz)

<0.0016 min (0.016s response time) (160 Hz)

>0.0016 min (0.031s response time) (160 Hz)

>0.0031 min (0.063 s response time) (80 Hz)

>0.0063 min (0.13 s response time) (40 Hz)

>0.013 min (0.25 s response time) (20 Hz)

>0.025 min (0.5 s response time) (10 Hz)

>0.05 min (1.0 s response time) (5 Hz)

>0.10 min (2.0 s response time) (2.5 Hz)

>0.20 min (4.0 s response time) (1.25 Hz)

>0.40 min (8.0 s response time) (0.62 Hz)

>0.85 min (16.0 s response time) (0.31 Hz)

Stoptime

As P

min

Advanced

Spectrum

Store: All

Range from: 190.0 to 400.0 nm

Step: 2.0 nm

Threshold: 10.0 mAU

Analog Output

Output 1:

Zero Offset: 5 %

Attenuation: 1000 mAU

Margin for negative Absorbance

100 mAU

Slit

4 nm

Autobalance

Prerun


Postrun

Lamps on required for acquisition

UV Lamp

# Detectors

For narrower peaks, increase data rates!  
Maintain Resolution at High Analysis Speed



## 80Hz versus 10Hz (20Hz) Data Rate

- Peak Width: - 55% (- 30%)
- Resolution: + 90% (+ 30%)
- Peak Capacity: + 120% (+ 40%)
- App. Column Eff.: + 260% (+ 70%)

Data Rate	Peak Width	Resolution	Peak Capacity
80 Hz	0.300	2.25	60
40 Hz	0.329	2.05	55
20 Hz	0.416	1.71	45
10 Hz	0.666	1.17	29
5 Hz	1.236	0.67	16


Sample: Phenones Test Mix  
Column: Zorbax SB-C18, 4.6x30, 1.8um  
Gradient: 50-100%ACN in 0.3min  
Flow Rate: 5ml/min

# Detector – Baseline noise

## Data rate

**Data rate = 10 Hz**


$$\text{Noise}_{(\text{PtoP}, 0.5 - 1.5 \text{ min})} = 1.48 \times 10^{-2}$$



File size<sub>(DAD1A.ch)</sub> = 9 KB


**Data rate = 160 Hz**

$$\text{Noise}_{(\text{PtoP}, 0.5 - 1.5 \text{ min})} = 5.43 \times 10^{-2}$$



File size<sub>(DAD1A.ch)</sub> = 45 KB

# Total Signal with Diode Array Detection



$$\frac{\text{Absorbance}_{\text{Sample wavelength}} + \text{Absorbance}_{\text{averaged over Bandwidth}}}{\# \text{ of wavelengths used}} + \frac{\text{Absorbance}_{\text{ref wavelength}} + \text{Absorbance}_{\text{averaged over Bandwidth}}}{\# \text{ of wavelengths used}} = \text{Total Absorbance}$$

# Care of Detector Flow Cells

Avoid the use of alkaline solutions with  $\text{pH} > 9.5$  which can attack quartz and impair optical performance.

Prevent crystallization of buffers or salts which will lead to blockage and damage.

Aqueous solvents can allow algae growth. Don't leave 100% water standing in the flow cell. When leaving LC idle, pump a mobile phase with at least 5-10% of organic solvent.

Observe the pressure limits of flow cells. Be careful when using detectors in series or fraction collectors.





# System Optimization

## Definitions of Critical Parameters


- Delay volumes
- Dispersion

# Typical HPLC System -Important Parameters





# Considerations for HPLC systems



## Gradient Delay or Dwell Volume

The volume between the point of mixing of solvents (usually in the mixing chamber or at the proportioning valves in the liquid chromatograph) and the head of an LC column.

## Extracolumn Volume

The volume between the effective injection point and the effective detection point, excluding the part of the column containing the stationary phase. It comprises the volumes of the injector, connecting lines and frits, and the detector. It determines the *extracolumn effects*.

# Conclusions

This Presentation was intended to be an introduction to The High Performance Liquid Chromatography Instrument. The HPLC is a widely used and powerful tool in analytical chemistry.

We discussed:

Typical HPLC Instrumentation

Understanding the Flow Path

- Pump
- Autosampler
- Column Compartment
- Common HPLC Detectors

General Maintenance and Housekeeping

General Settings

Resources for Support

# Resources — Primers

[5990-7595EN](#)

## The LC Handbook

Guide to LC Columns and Method Development

[5991-2359EN](#)

## Two Dimensional Liquid Chromatography

[5990-3777EN](#)

## High Performance Capillary Electrophoresis

[5991-5509EN](#)

## Supercritical Fluid Chromatography

[5989-6639EN](#)

## Principles in Preparative HPLC

[5991-3326EN](#)

## Sample Preparation Fundamentals for Chromatography


[5980-1397EN](#)

## Fundamentals of UV-visible Spectroscopy



# Resources for Support

- Collection of LC resources: [https://community.agilent.com/docs/DOC-1852-lc-insights-to-go#jive\\_content\\_id\\_LC\\_Troubleshooting](https://community.agilent.com/docs/DOC-1852-lc-insights-to-go#jive_content_id_LC_Troubleshooting)
- Agilent support resources: <https://community.agilent.com/community/resources>
- Agilent University: <http://www.agilent.com/crosslab/university>
- Agilent resource center: <http://www.agilent.com/chem/agilentresources>
- InfinityLab Supplies Catalog ([5991-8031EN](tel:5991-8031EN))
- Your local FSE and Specialists
- Youtube – [Agilent Channel](#)
  
- Sales and support phone assistance (US and Canada):  
1-800-227-9770 [Phone Tree Navigation Assistance](#)



[gc-column-support@agilent.com](mailto:gc-column-support@agilent.com)  
[lc-column-support@agilent.com](mailto:lc-column-support@agilent.com)  
[spp-support@agilent.com](mailto:spp-support@agilent.com)  
[spectro-supplies-support@agilent.com](mailto:spectro-supplies-support@agilent.com)





# Thank you!!!!



# Questions?