

Here's Proof: This Ultra Inert GC Column Outperforms the Competition


Agilent J&W DB-5ms Ultra Inert GC columns



New study confirms consistent peak shape for active compounds—even after extended periods of time at maximum temperature

The Agilent J&W Ultra Inert GC column family pushes industry standards for consistent column inertness and exceptionally low column bleed. The result? Lower detection limits and more accurate data for active analytes.

In this study, we explored how a sustained maximum programming temperature of 350 °C would affect the peak shape of 1,2-butanediol. The peak shape symmetry stayed near 1.0 for over 50 hours, demonstrating the robustness and superior inertness of DB-5ms UI columns.



Peak	Name
1	1,2-Butanediol
2	Decane
3	1-Octanol
4	2,6-Dimethylphenol
5	2,6-Dimethylaniline
6	Naphthalene
7	1-Decanol
8	Tridecane
9	Methyl Decanoate

Conditions

GC System: Agilent 8890/FID
Carrier: Helium, constant flow, 3.3 mL/min
Column: J&W DB-5ms Ultra Inert
30m x 0.25mm x 0.25 μ m
(122-5532UI)
Oven: 120 °C
FID: 260 °C, hydrogen: 40 mL/min; air:
400 mL/min
Inlet: Split 200:1; 250 °C
Injection: 1 μ L

To learn how Agilent Inert Flow Path solutions can give you the utmost confidence in your results, visit

www.agilent.com/chem/inert

DE44389.682777778

This information is subject to change without notice.

© Agilent Technologies, Inc. 2021
Published in the USA, July 19, 2021
5994-3757EN