


# GC/FID & GC/MS RTL Flavor Databases

**User Contributed by**  
**Agilent Technologies**  
**Wilmington, Delaware, USA**




**Agilent Technologies**  
Innovating the HP Way

# Why Use Retention Time Locking ?

---


Retention times shift with column maintenance or detector



# Using Retention Time Locking

---

Retention times are reproducible to 0.030 min or better



# Benefits of RTL

---

- Pattern profiles can be superimposed
- Quantitation database times stay the same
- SIM group switching times stay the same
- Integration event times stay the same
- Troubleshooting simplified
- Comparing data from other testing facilities
- GC/FID database can search on retention time instead of retention index
- GC/MS Screener database can search on retention time as a qualifier

# Database Development

---

- **Developed by**
  - **Research Institute for Chromatography, Kortrijk, Belgium**
  - **Quest International, Naarden, The Netherlands**
  - **Laboratory of Organic Chemistry, University of Gent, Belgium**
  - **In cooperation with Agilent Technologies**
- **All retention time and spectral data collected on Agilent GC/FID and GC/MSD systems**

# There are 2 GC/MS Databases

---

- **Retention Time Locked Screener Database**
  - Flavor2.scd
  - 409 compounds with precise retention times
  - target ion and up to 3 qualifier ions
- **Mass Spectral Library**
  - Flavor2.L
  - GC/MS Spectra for the 409 compounds
- **RTLlocked to n-pentadecane at 27.500 min**
  - 30 m x 250  $\mu\text{m}$  x 0.25  $\mu\text{m}$  HP-5MS, 19091S-433

# **There is also 1 GC/FID Database**

---

- **Retention Time Locked Database**
  - **Flavfid.rtt**
  - **Same 409 compounds as the GC/MS databases**
  - **precise retention times**
- **RTLlocked to n-pentadecane at 70.000 min**
  - **50 m x 320  $\mu$ m x 1.05  $\mu$ m HP-5, 19091J-215**


# Using the GC/FID Flavor Database

---

- Flavfid.rtt is copied to the C:\Hpchem\Rtl folder
- The GC/FID system is retention time locked and an RTL method is established
- Samples are run using the RTL method
- Sample datafile is analyzed using Flavfid.rtt
  - interactively with the RTSearch menu item
  - automatically with the user contributed macro Autosrch.mac
- Results can be saved or printed



# Spanish Orange Oil - RTLocked GC/FID



# Search Results - Spanish Orange Oil

---

Peak number	GC-FID <u>t<sub>R</sub></u> (min)	t <sub>R</sub> identification	GC-MS t <sub>R</sub> (min)	MS + t <sub>R</sub> identification
1	<u>26.793</u>	$\alpha$ -pinene	5.172	$\alpha$ -pinene
2	30.042	1-octen-3-ol 3-(methylthio)-1-propanol sabinene	6.181	sabinene
3	30.539	hexanoic acid $\beta$ -pinene 6-methyl-5-hepten-2-one	6.282	$\beta$ -pinene
4	31.053	2-octanone myrcene furfuryl acetate	6.658	myrcene
5	31.987	octanal	6.987	octanal
6	33.190	trans 2-hexenoic acid $\Delta$ -3-carene	7.267	$\Delta$ -3-carene
7	35.001	limonene benzylalcohol ocimene	8.130	limonene

# Search Results - Spanish Orange Oil

---


Peak number	GC-FID $t_R$ (min)	$t_R$ identification	GC-MS $t_R$ (min)	MS + $t_R$ identification
8	40.162	n-undecane cis-3-hexenylpropionate $\delta$ -hexalactone 1-methyl-2,3-cyclohexadione linalool methyl benzoate	10.391	linalool
9	48.728	dihydrocarveol methyl salicylate estragole n.decanal octylacetate	14.750	n-decanal
10	71.366	anisylpropionate valencene piperonyl acetate	27.134	valencene

# Using the GC/MSD Flavor Databases

---

- Flavor2.scd and Flavor2.L are copied to the Database folder, usually C:\Database
- The GC/MS system is retention time locked and an RTL method is established
- Samples are run using the RTL method
- Sample data file is Screened using Flavor2.scd
  - **Screeener report is generated**
- **Screeener Results are reviewed/modified**
- **Peak spectra can be searched/compared using Flavor2.L**

# Spanish Orange Oil - RTLocked GC/MSD



# Search Results - Spanish Orange Oil

---

Peak number	GC-FID $t_R$ (min)	$t_R$ identification	GC-MS $t_R$ (min)	MS + $t_R$ identification
1	26.793	$\alpha$ -pinene	5.172	$\alpha$ -pinene
2	30.042	1-octen-3-ol 3-(methylthio)-1-propanol sabinene	6.181	sabinene
3	30.539	hexanoic acid $\beta$ -pinene 6-methyl-5-hepten-2-one	6.282	$\beta$ -pinene
4	31.053	2-octanone myrcene furfuryl acetate	6.658	myrcene
5	31.987	octanal	6.987	octanal
6	33.190	trans 2-hexenoic acid $\Delta$ -3-carene	7.267	$\Delta$ -3-carene
7	35.001	limonene benzylalcohol ocimene	8.130	limonene

# Search Results - Spanish Orange Oil

---

Peak number	GC-FID $t_R$ (min)	$t_R$ identification	GC-MS $t_R$ (min)	MS + $t_R$ identification
8	40.162	n-undecane cis-3-hexenylpropionate $\delta$ -hexalactone 1-methyl-2,3-cyclohexadione linalool methyl benzoate	10.391	linalool
9	48.728	dihydrocarveol methyl salicylate estragole n.decanal octylacetate	14.750	n-decanal
10	71.366	anisylpropionate valencene piperonyl acetate	27.134	valencene

# Screenener Variables

---

- **The number of probable and possible Hits can be globally changed through these variables**
  - **Extraction Window for retention time**
  - **Qualifier Mode, relative or absolute**
  - **Qualifier %**
  - **Zero qualifiers, included or excluded**
  - **Subtraction Mode**
- **Integration parameters define peak detection**



# Benefits of Using RTL Databases

---

- **Fastest identification compounds**
  - Screen for 409 compounds in < 15 sec
- **Fastest confirmation of compounds**
  - Eliminate hits with wrong r.t.
- **Precise reproducible r.t. on GC and GC/MS**
  - No need to update calibration tables
- **Compounds identified by both retention time and spectral information**
- **No additional cost for user-contributed databases**

# Flavor Database - RT Order 1 of 6

---

Ethyl alcohol	n-Propyl acetate	4-Hexene-3-one
Isopropyl formate	Ammonium sulfide	5-Methyl-2,3-hexanedione
Isopropyl alcohol	Methyl butyrate	Ethylcrotonate; ethyl-trans-2
Isobutyraldehyde	Isoamyl alcohol	Ethyl-2methylbutyrate
Diacetyl	4-Methyl-2-pentanone	Ethyl isovalerate
Butyraldehyde; butanal	Isobutyric acid	trans-2-Hexenal
Methyl ethyl ketone	Ethyl isobutyrate	Furfuryl alcohol
Hexane	Amyl alcohol	cis-3-Hexenol
Ethyl acetate	Isobutylacetate	2-Hexen-1-ol (trans)
isobutyl alcohol	Butyric acid	1-Hexanol
3-Methylbutyraldehyde; iso	Methyl-2-methylbutyrate	Isobutylpropionate
Butylalcohol	Isoamyl formate	Isoamyl acetate
Isopropyl acetate	3,4-Hexanedione	2-Methylbutylacetate
2-Methylbutyraldehyde; 2-	n-Octane	Valeric acid
Propionic acid	Ethyl butyrate	Styrene
1-Penten-3-ol	4-Methyl-3-penten-2-one	Isopropyl isovalerate
Methyl isobutyrate	Hexenal	n-Nonane
2-Pentanone	2-Methyltetrahydrofuran-3-one	2-Heptanol
Isobutyl formate	Propyl propionate	cis-4-Heptenal
Acetyl propionyl	n-Butyl acetate	Ethyl valerate
2-Pentanol	Ethyl lactate	Heptanal
n-Valeraldehyde	Methylvalerate	Butyl propionate
Pyruvic acid; 2-oxopropan	2-Methylpyrazine	trans,trans-2,4-Hexadienal
Heptane	2-Methyl-2-pentenal	2-Furyl-methylketone; 2-acetyl
Ethyl propionate	Isovaleric acid	gamma-Butyrolactone

# Flavor Database - RT Order 2 of 6

---

Isobutyl isobutyrate	2-Octanone	p-Cymene
Anisole	Myrcene	2-Ethyl-1-hexanol
Methyl caproate; methyl h	3-Octanol	Limonene
Hexyl formate	Furfuryl acetate	Eucalyptol
alpha-Pinene	Butyl butyrate	Benzyl alcohol
Ethylacetoacetate	n-Decane	2-Isobutylthiazole
4-Methylpentanoic acid	3-Hexenoic acid	Isopropyl hexanoate
Propyl isovalerate	2-Octanol	3,7-Dimethyl-1,3,6-octatriene
Camphene	2,3,5-Trimethylpyrazine	Cyclohexyl acetate
gamma-Valerolacton	Octanal	n-Butyl-2-methylbutyrate
Isobutyl butyrate	alpha-phellandrene	Phenylacetaldehyde
trans-2-Heptenal	Isobutyl isovalerate	Ethyl-trans-2-hexenoate
Benzaldehyde	cis-3-Hexenyl acetate	trans-2-Hexenoic acid
Heptyl alcohol	Ethyl-3-hexenoate	Butyl isovalerate
Isoamyl propionate	trans,trans-2,4-Heptadienal	o-Cresol
2-Methoxy-3-methylpyrazin	delta-3-carene	gamma-Hexalactone
2-Ethylbutyl acetate	Hexyl acetate	2,6-Dimethyl-5-heptanal
sabinene	trans-2-Hexen-1-yl-acetate	3,4-Dimethyl-1,2-cyclopentadienol
Methyl-2-furoate	2-Methyl-2-pentanoic acid	Isoamylbutyrate
3-(Methylthio)-1-propanol	1,4-Cineole (Isocineole)	Amyl butyrate
1-octen-3-ol	trans-2-Hexen-1-yl acetate	trans-2-octenal
Phenol	alpha-terpinene	gamma-Terpinene
beta-Pinene	p-Methylanisole	alpha-Methylbenzylalcohol
Hexanoic acid	2-Acetylpyrazine	Levulinic acid
6-Methyl-5-hepten-2-one	Methyl-3-(methylthio)propionate	Ethyl levulinate

# Flavor Database - RT Order 3 of 6

---

Methyl-2-pyrrolyl ketone; 2	Nonanal	Menthone
Acetophenone	6-Methyl-3,5-heptadien-2-one	Benzyl acetate
Diethyl malonate	3-Acetylpyridine	Isoborneol (Isomer 2)
p-Cresol	Isoamyl isovalerate	Benzoic acid
Linaloloxide (cis, isomer B	Methylbutyl-2- isovalerate	p-Methylphenyl acetate
m-Cresol	Rose oxide (cis); tetrahydro-4-m	Ethylbenzoate
Benzyl formate	Phenethyl alcohol	cis-6-Nonen-1-ol
Heptanoic acid	Heptyl acetate	Nonanol
Allyl hexanoate	Fenchyl alcohol	L-Menthol
Methyl phenylethyl ether	Isophorone	Octanoic acid; caprylic acid
3,5-Dimethyl-1,2-cyclopent	Ethyl-3-hydroxyhexenoate	Phenethyl formate
2,3,5,6-Tetramethylpyrazin	Methylnicotinate	4-Carvomenthenol; terpinen
Terpinolene	2-Methylheptanoic acid	Methyl phenyl acetate
Linaloxide (trans, isomer A	Homofuronol; 2-ethyl-4-hydroxy	Allyl heptanoate
Guaiacol	Isopulegol	Diethyl succinate
2-Nonanone	Isobutyl caproate; isobutyl hexa	4'-Methylacetophenone
1-Methyl-2,3-cyclohexadio	Citronellal	2,4-Dimethylbenzaldehyde
delta-Hexalactone	gamma-Heptalactone	cis-3-Hexenyl butyrate
Methyl benzoate	trans-2,cis-6-Nonadienal	alpha-Terpineol
Ethyl heptanoate	Neroloxide	Butyl capronate; butyl hexar
n-Undecane	Isoborneol (isomer 1)	Hexyl butyrate
Linalol	alpha,alpha-Dimethylphenethyl a	Methylsalicylate
Phenethylamine	trans-2-Nonenal	Dihydrocarveol
cis-3-Hexenylpropionate	3-Phenylpropionaldehyde	cis-4-Decenal
cis-6-Nonenal	2-Methyl-4-propyl-1,3-oxathiane	Propyl heptanoate

# Flavor Database - RT Order 4 of 6

---

alpha-Methylbenzylacetate	Isoamyl hexanoate	Carvacrol
3,7-Dimethyl-1-octanol	p-Anisaldehyde; p-methoxybenz	Methyl octine carbonate
Ethyl octanoate	d-Piperitone	Geranyl formate
Estragole	Geraniol	5-Methylquinoxaline
Safranal	Phenethyl acetate	Undecanal
n-Dodecane	Phenylacetic acid	Nonyl acetate
cis-3-Hexenyl lactate	gamma-Octalactone	trans,trans-2,4-Decadienal
n-decanal	Benzyl propionate	Citronellic acid; 3,7-dimethy
Octyl acetate	Cinnamaldehyde	alpha,alpha-Dimethyl phenol
trans,trans-2,4-Nonadiena	Ethyl salicylate	Isobutyl benzoate
Linalyl formate	Citral	Dihydrocarvylacetate
Phenylacetaldehyde dimet	1-Decanol	Piperonal; heliotropin
beta-Cyclocitral	Isopropyl phenylacetate	Anisyl formate
Cyclohexyl butyrate	Nonanoic acid	Carvyl acetate
Cuminaldehyde	Anisyl alcohol	Propylphenyl acetate
Nerol, cis-geraniol	delta-Octalactone	Methyl anthranilate
Citronellol	trans-Anethole	3-Phenylpropionic acid
3-Phenyl-1-propanol	Butyl heptanoate	Benzyl butyrate
cis-3-Hexenyl-2-methyl-bu	2-undecanone	Terpinyl acetate
Quinoline	Heptyl butyrate	2,6-Dimethoxyphenol
cis-3-Hexenyl isovalerate	Indole	Acetanisole
4-Methyloctanoic acid	Thymol	Phenethyl propionate
Hexyl-trans-2-butenate	Menthyl acetate; d,L-methyl-2-(m	Triacetin
Carvone	Ethyl Nonanoate	Benzylidene acetone
Ethylphenyl acetate	Benzyl isobutyrate	Eugenol

# Flavor Database - RT Order 5 of 6

---

Butylbutyryllactate	beta-Caryophyllene	beta-Ionone
gamma-Nonalactone	Anisyl acetate	Veratraldehyde; methyl vanil
Neryl acetate	Linalyl butyrate	Butylated hydroxyanisole
3-Phenylpropyl acetate	Isobutyl salicylate	Dimethyl benzyl carbonyl but
Decanoic acid	Allyl cyclohexyl propionate	Valencene
Dihydrocoumarin	trans-Cinnamic acid	delta-Decalactone
Methyl cinnamate	Geranylacetone (Isomer 1); (E)-6,1	Pentadecane
Dihydrojasmane	Ethyl methylphenylglycidate (iso	Piperonyl acetate
Damascenone	Octyl isovalerate	Anisyl propionate
Geranyl acetate	Phenethyl butyrate	alpha-Irone (isomer 1) (E-); r
Hexyl hexanoate	Citronellyl propionate	Ethyl-3-phenylglycidate
delta-Nonalactone	Cinnamyl acetate	Methyl-gamma-ionone (isom
Octyl butyrate	Ethyl-p-anisate	Propenyl guaethol
Benzyl isovalerate	Isoeugenol	Ethyl methylphenylglycidate
Phenethyl isobutyrate	Geranylacetone(isomer 2); (E)-6,1	Citronellyl butyrate
Vanillin	6-Pentyl-alpha-pyrone; 5-hydrox	Nerolidol, (Z-); nerolidol, (cis
Ethyl laurate; ethyl decano	Ethyl cinnamate	alpha-Irone (isomer 2); meth
cis-Jasmone	gamma-Decalactone	nonyl isovalerate
Diphenylether	Ethyl-trans-2, cis-4-decadienoate	Bisabolene
n-Tetradecane	Lauryl alcohol	6-Methylcoumarin
Dimethyl anthranilate	Geranyl propionate	Geranyl butyrate
Lauric aldehyde	5-Hydroxy-2-decanoic acid delta	Nerolidol, (E-); nerolidol, (tra
Decyl acetate	Methyl-gamma-Ionone (isomer 1	Lauric acid
Ethyl anthranilate	Isoamyl octanoate	cis-3-Hexenylbenzoate
beta-Damascone	Isoamyl octanoate	Isoamyl salicylate

# Flavor Database - RT Order 6 of 6

---

gamma-Undecalactone

Cinnamyl isobutyrate

Ethyl laurate

n-Hexadecane

Anisyl butyrate

delta-Undecalactone

Lauryl acetate

Methyljasmonate

Methyl dihydrojasmonate

gamma-Dodecalactone

Farnesol (Z,E-)

delta-Dodecalactone

Farnesol (E,E-)

Amberoxid; dodecahydro-

Benzyl benzoate

Myristic acid

Ethyl myristate

Guaiyl acetate

n-Octadecane

Nootkatone

Isopropyl myristinate

omega-Pentadecalactone

Phenethyloctanoate

Phenethyl phenyl acetate

omega-6-Hexadecalactone

Palmitic acid

Ethyl palmitate

n-Eicosane

Benzyl cinnamate

Ethyl oleate

Stearic acid

ethyloctadecanoate, Ethylsteara

Methional; 3-(methylthio)propion

n-Tetracosane

# Agilent Instrumentation Required, GC/MSD

---

- **Gas Chromatograph**

- **6890N with autosampler (tray & injector) and split/splitless inlet**

- **Mass Spectrometer**

- **5973N with EI source**

- **Software**

- **GC/MSD ChemStation G1701CA version C.00.01 or higher (includes RTL and Screener)**

- **Column**

- **HP-5MS 30m x 250 $\mu$ m x 0.25 $\mu$ m part # 19091S-433**



# Agilent Instrumentation Required, GC/FID

---

- **Gas Chromatograph**

- **6890N with autosampler (tray & injector), FID and split/splitless inlet**

- **Software**

- **GC ChemStation revision A.08.03 or higher**
- **GC RTL revision A.06.01 or higher**

- **Column**

- **HP-5 50m x 320 $\mu$ m x 1.05 $\mu$ m part # 19091J-215**

# GC/FID Flavor Parameters

---

**GC** Agilent 6890N

**Oven**

Ramp	'C/min	'C	Hold min
------	--------	----	----------

Initial		50	0.00
---------	--	----	------

Ramp 1	2	280	0.00
--------	---	-----	------

**Runtime** 115 min

**Inlet** Split/Splitless

**Temp** 250 °C

**Mode** Split

**Pressure** 12.33 psi (adj to lock)

**Split Ratio** 25:1 for 5% solution  
100:1 for oils

**Column** HP- 5 part # 19091J-215

**Mode** Constant pressure

**Pressure** 12.33 psi

**Initial Flow** 1.6 mL/min

**Outlet Pressure** Ambient

**Detector**

**FID**

**Temp** 250 °C

**Hydrogen** 40 mL/min

**Air** 450 mL/min

**Nitrogen** 45 mL/min

**Makeup**

**Injector**

**Sample Washes** 0

**Sample Pumps** 3

**Injection volume** 1 uL for 5 % soln.  
0.2 uL for oils

**Solv A/B washes** 2 each

**Viscosity Delay** 0

**Plunger Speed** Fast

# GC/MSD Flavor Parameters

---

**GC** Agilent 6890N

**Oven**

Ramp	'C/min	'C	Hold min
------	--------	----	----------

Initial		60	0.00
---------	--	----	------

Ramp 1	3	240	0.00
--------	---	-----	------

Runtime	60 min		
---------	--------	--	--

**Inlet** Split/Splitless

**Temp** 250 'C

**Mode** Split

**Pressure** 9.43 psi (adj to lock)

**Split Ratio** 25:1 for 5% solution  
100:1 for oils

**MSD** Agilent 5973N

**Solvent Delay** 2.0 min

**Scan Range** 40 to 400

**Threshold** 150

**Sampling** 2

**Quad Temp** 150 'C

**Source Temp** 230 'C

**Transfer Line** 280 'C

**Column** HP- 5MS part # 19091S-433

**Mode** Constant pressure

**Pressure** 9.43 psi

**Initial Flow** 1.2 mL/min

**Outlet** MSD

**Outlet Pressure** Vacuum

**Injector**

**Sample Washes** 0

**Sample Pumps** 3

**Injection volume** 1 uL for 5 % soln.

0.2 uL for oils

**Solv A/B washes** 2 each

**Viscosity Delay** 0

**Plunger Speed** Fast

# Recommended Reading

---

- “Analysis of Essential Oil Compounds Using Retention Time Locked Methods and Retention Time Databases”  
publication# 5988-6530EN
- Direct your web browser to [www.agilent.com/chem](http://www.agilent.com/chem)
- Click on the **Library** link in the left navigation panel
- In the literature search section, enter **5988-6530EN** as the keyword and click the **Search** button

# More Recommended Reading

---

- Retention Time Locking with the HP G1701BA MSD Productivity ChemStation (5968-3433)
- Efficient Screening for Pesticides and Endocrine Disrupters Using the HP 6890/ 5973 GC/MSD System (5968-4884)
- Retention Time Locking: Concepts and Applications (5966-2469)
- Precise Time-Scaling of Gas Chromatographic Methods Using Method Translation and Retention Time Locking (5967-5820)
- Direct your web browser to [www.agilent.com/chem](http://www.agilent.com/chem)
- Click on the **Library** link in the left navigation panel
- In the **keyword** field, type in the 8 digit number as it is listed above for the publication you would like to view