

# The Agilent 7696A Sample Prep WorkBench

Techniques and Strategies to  
Build Automated Sample Prep  
Methods

Jim McCurry  
Senior Scientist  
Little Falls Site  
Wilmington, Delaware USA

# Outline

- What is the 7696A WorkBench
- Steps to translating and building WorkBench methods from manual methods
- Optimizing WorkBench method conditions
- Results using WorkBench Methods

# Agilent 7696A Features and Capabilities

## Simple Sample Prep:

- **Dilution** / Aliquoting / Reconstitution
- **Additions** (standards, reagents, etc)
- **Heating** (Reactions, Derivatization, etc)
- **Mixing** – Vortex


## Additional Capabilities:

- Sample Tray Heating – 50 vial positions\*
- Sample Tray Peltier Cooling – 50 vial positions\*
- Dedicated Sample Prep Software
- Liquid/Liquid Extraction




\* Option / Accessory

# Typical Manual Sample Prep – ASTM D6584



# WorkBench Sample Prep – ASTM D6584



# Translating Manual Methods to WorkBench

## IP 585 – GC/MS Analysis of Trace FAME in Jet Fuel

- Method measures 5 to 100 mg/kg total FAME contamination in Jet Fuel
- GC/MS analysis technique requires and expensive, deuterated internal standard
  - Methyl Heptadecanoate-d33 (C17:0d33)
- Step 1 – Convert Manual method to a step-by-step script
- Step 2 – Scale each step to fit into 2-mL ALS vials
- Step 3 – List resources need for each step
- Step 4 – Layout resources in WorkBench software
  - Chemicals, standards, solvents, reagents, syringes, wash solvents, waste vials
- Step 5 – Build method using WorkBench software
- Step 6 – Optimize WorkBench conditions
  - Syringe conditions, mixing, heating or cooling

# Translating Manual Methods to WorkBench

## Step 1 – Convert Manual method to a step-by-step script

Convert manual method to a “Script” of steps

### Steps Manual Prep of 10 mg/kg IP585 Standard

- 1 Add 900 uL dodecane to vial
- 2 Add 100 uL working standard solution (100 mg/kg)
- 3 Add 10 uL ISTD solution (1000 mg/kg C17:0d33)
- 4 Mix

### Steps Manual Prep of Sample

- 1 Add 1000 uL of sample to vial
- 2 Add 10 uL ISTD solution (1000 mg/kg C17:0d33)
- 3 Mix

FAME IN AVTUR – GC-MS- IP 585

6.3.1 Used containers are permitted provided it can be confirmed they have only been used for aviation turbine fuel containing < 5 mg/kg FAME.

NOTE - New sample containers are strongly recommended due to concerns over the difficulty in removing all traces of FAME retained from previous samples.

6.3.2 Rinse all sample containers and their closures with the product to be sampled at least three times. Each rinse shall use product with a volume of 10% to 20% of the container volume. Each rinse shall include closing and shaking the container for a minimum of 5 seconds and then draining the product.

7. Internal Standard and Calibration Standards preparation

7.1 Internal Standard Solution (ISS) preparation: Prepare a stock internal standard solution containing nominally 1000 mg/L of methyl heptadecanoate-d33 internal standard as follows.

7.1.1 Place a 100 ml volumetric flask (5.11) onto the balance (5.3) and tare. Using a dropping pipette (5.6), add 100mg +/- 0.5mg of methyl heptadecanoate-d33 (4.5) and make up to 100ml with dodecane (4.3) to give a 1000mg/L solution.

7.1.2 Store 1000 mg/L standard in a tightly closed glass container in a refrigerator held at 4° C +/- 2° C when not in use and use within three months. Before use, examine for any phase separation or discoloration; shake vigorously, and let stand to allow for removal of air bubbles. Discard the standard if it shows sediment, phase separation or discoloration.

7.2 Bulk Calibration Solution (BCS) preparation: Prepare a bulk calibration solution (BCS) containing nominally 1000 mg/kg each of all of the FAME reference compounds listed in Table 1 in dodecane.

7.2.1 Place the conical flask (5.4) onto the balance (5.3) and tare. Using a dropping pipette (5.6) place 0,1000g ± 0,0010 g of each FAME species into the flask and record the masses added. Add dodecane (4.3) to give a total mass of 100 g ± 0,05 g; Stopper the flask and mix the contents.

7.2.2 Transfer the BCS to a tightly closed glass container and store in a refrigerator held at 4° C +/- 2° C when not in use and use within three months. Before use, examine for any phase separation or discoloration; shake vigorously, and let stand to allow for removal of air bubbles. Discard the standard if it shows sediment, phase separation or discoloration.

7.2.3 Calculate the stock solution concentration for each FAME component, (Y) using equation (1).

$$\text{Conc (mg/kg) of FAME species Y in BCS} = \frac{\text{mass of FAME species Y (g)} \times 10^4}{\text{total mass of BCS (g)}} \quad (1)$$

7.3 Working Standard Solution (WSS) and Working Calibration Standards

7.3.1 Prepare a nominal 100 mg/kg working standard solution (WSS) by pipetting 1000 µl of BCS into the 10ml volumetric flask and dilute to 10 ml with dodecane. Calculate the working standard solution (WSS) concentration for each FAME component (Y) using equation (2).

$$\text{Conc (mg/kg) of FAME species Y in WSS} = \frac{\text{conc of FAME species Y in BCS}}{10} \quad (2)$$

7.3.2 Prepare a set of working calibration standards containing nominally 2, 4, 6, 8, 10, 20, 40, 60, 80 and 100 mg/kg of each reference compound by dilution of the WSS (7.3.1).

Table 2: Volumetric Dilutions for preparing standards in the range 0 – 100mg/kg from the 100mg/kg WCS

Nominal Std concentration (mg/kg)	Volume WSS (µl)	Volume Dodecane (µl)	Volume Internal Standard (µl)
100	1000	0	10
80	800	200	10
60	600	400	10
40	400	600	10
20	200	800	10
10	100	900	10
6	60	940	10
4	40	960	10
2	20	980	10
0	0	1000	10

Issued under licence to ConocoPhillips employees only.

IMPORTANT: This file is subject to a licence agreement issued by the Energy Institute, London, UK. All rights reserved. It may only be used in accordance with the licence terms and conditions. It must not be forwarded to, or stored or accessed by any unauthorised user. Enquiries: e: pubs@energyinst.org.uk; +44 (0)207 467 7100



# Scaling Manual Prep for 2-mL WorkBench vials

## Step 2 – Scale each step to fit into 2-mL ALS vials

Steps	Manual Prep	10:1 Scaling -->	WorkBench Prep
1	Add 900 uL of dodecane to a vial		Add 90 uL of dodecane to a vial
2	Add 100 uL of working calibration standard		Add 10 uL of working calibration standard
3	Add 10 uL ISTD solution		Add 1 ISTD solution
4	Mix		Mix

- Manual sample preps must be scaled to fit 2-mL vial
- Reduce amounts of solvents and expensive chemicals
- High recovery vials needed to contain 100 uL of final solution volume





# Translating Manual Methods to WorkBench

## Step 3 – List resources need for each step

- Four chemical resources are needed for this sample prep
- Syringe size, draw speeds, and dispense speed are optimized for each type of chemical resource

<b>Resource</b>	<b>Tray Position</b>	<b>Syringe Size</b>	<b>Syringe Draw Speed</b>	<b>Syringe Dispense Speed</b>
Dodecane	51 - 52	250 uL	2500 uL/min	5000 uL/min
Stock Standard (100 mg/kg each FAME)	71 -72	25 uL	250 uL/min	500 uL/min
ISTD (1000 mg/kg C17:0d33)	91	25 uL	250 uL/min	500 uL/min

# Translating Manual Methods to WorkBench

## Step 4 – Layout resources in WorkBench software

Sample Prep Resource Layout Editor Version 3.1.34.0

Tray Resources Wash/Waste Vial Assignment

Resource Name: Internal Standard  
Resource Type: Chemical Resource

Use Type:  By Volume  
Useable Volume per Vial (µL): 500  
 By Use  
Uses per Vial: 1

Display Color: Maroon

Default Syringe Parameters

For Syringe Size (µL): 25  
Wash Volume (µL): 5  
Pump Volume (µL): 20  
Draw Speed (µL/min): 250  
Dispense Speed (µL/min): 500  
Draw Needle Depth Offset (mm): -2.0  
Use Needle Depth Offset for Dispense:   
Viscosity Delay (s): 2  
Air Gap (% Syringe Size): 0  
Overfill (% Syringe Size): 5

Layout Comment:

Vial Range: 91

Add Remove Replace Cancel

Color	Name	Resource Type	Vial Range	Usage
Green	Dodecane	Chemical Resource	51-52	1500 µL
Red	100 ppm Stock Std	Chemical Resource	71-72	1500 µL
Blue	Internal Standard	Chemical Resource	91	500 µL

Save Layout Close

Sample Prep Resource Layout Editor Version 3.1.34.0

Tray Resources Wash/Waste Vial Assignment

Front Turret

Solvent A  
Solvent A Vials: A1-A6

Solvent B  
Solvent B Vials: B1-B3

Default Syringe Parameters

For Syringe Size (µL): 250  
Wash Volume (µL): 100  
Draw Speed (µL/min): 625  
Dispense Speed (µL/min): 1250  
Viscosity Delay (s): 0

Waste A Vials: WA1-WA3 Waste B Vials: WB1-WB2

Back Turret

Solvent A  
Solvent A Vials: A1-A6

Solvent B  
Solvent B Vials: B1-B3

Default Syringe Parameters

For Syringe Size (µL): 25  
Wash Volume (µL): 10  
Draw Speed (µL/min): 750  
Dispense Speed (µL/min): 1500  
Viscosity Delay (s): 2



Waste A Vials: WA1-WA3 Waste B Vials: WB1-WB2

Save Layout Close

# Translating Manual Methods to WorkBench

## Step 5 – Build method using WorkBench software

- WorkBench software builds methods graphically
- “Drag and drop” WorkBench functions into the method program



# Translating Manual Methods to WorkBench

## Step 5 – Build method using WorkBench software

Setup Method

Agilent 7696A Sample Prep Method | Agilent 7696A Configuration

Import Export

Process in Batch Mode

Version 3.1.34.0

Actions

Add Mix Heat Wait Flag as result Move Wash

Steps

1. Wash with 200  $\mu\text{L}$  of Front Solvent A 1 times at Front Tower  
2. Add 90  $\mu\text{L}$  of Dodecane to Sample at Front Tower (washes, pumps)

Program

1. Wash → 2. Add

**Edit Add Step**

Add: 90  $\mu\text{L}$

Of: Dodecane

To: Sample

Check box to enable washes or pumps, then select repetitions.  
"Tower Selection" and "Dispense Settings" are always enabled.

Reset Selection to Defaults

Tower Selection  
 Solvent Prewash 1  
 Solvent Prewash 2  
 Dispense Wash  
 Dispense Pumps  
 Dispense Settings  
 Solvent Postwash 1  
 Solvent Postwash 2

**Dispense Settings**

Draw Speed ( $\mu\text{L}/\text{min}$ ): 1250

Dispense Speed ( $\mu\text{L}/\text{min}$ ): 5000

Draw Needle Depth Offset (mm): -2.0

Use Needle Depth Offset for Dispense:

Viscosity Delay (s): 2

Air Gap (% Syringe Size): 0

Overfill (% Syringe Size): 5

OK Cancel << Advanced

# Translating Manual Methods to WorkBench

## Step 5 – Build method using WorkBench software

**Steps**    **WorkBench Prep of 10 mg/kg IP585 Standard**

1	Solvent wash 250 µL syringe (front tower)
2	Add 90 µL dodecane to vial
3	Solvent wash 25 µL syringe
4	Add 10 µL working standard solution (100 mg/kg)
5	Solvent wash 25 µL syringe (back tower)
6	Add 1 µL ISTD solution (C17:0d33)
7	Mix

Resource Name	Resource Type	Usable Volume/Vi
ISTD	Chemical Resource	500 µL
Dodecane	Chemical Resource	1500 µL
100 ppm Stock Standard	Chemical Resource	1500 µL
Front Solvent A	Turret Location	2000 µL
Front Solvent B	Turret Location	2000 µL

# Translating Manual Methods to WorkBench

## Syringe Bubble Formation

### Step 6 - Optimize WorkBench conditions



This Guy Likes Bubbles




### WorkBench Does Not Like Bubbles

- Bubble formation in the syringe results in poor precision and accuracy when dispensing liquids

# What Causes Bubble Formation

- Cavitation
  - rapid change in pressure at the face of the plunger
  - low-boiling liquids
- Solvent degassing
  - i.e. dissolved air in methanol
- In-Vial Vacuum Formation
  - usually found when withdrawing large volumes from vial
  - overfilled vials
  - viscous liquids



# Techniques to Prevent Syringe Bubbles

Check box to enable washes or pumps, then select repetitions.  
"Tower Selection" and "Dispense Settings" are always enabled.

Reset Selection to Defaults

- Tower Selection
- Solvent Prewash 1
- Solvent Prewash 2
- Dispense Wash
- Dispense Pumps
- Dispense Settings
- Solvent Postwash 1
- Solvent Postwash 2

**Dispense Pumps**

Number of Pumps: 3

Pump Volume (µL): 400

Draw Speed (µL/min): 1250

Dispense Speed (µL/min): 2500

Needle Depth Offset (mm): 0

Viscosity Delay (s): 2

Use Dispense Pumps to blow out bubbles

Check box to enable washes or pumps, then select repetitions.  
"Tower Selection" and "Dispense Settings" are always enabled.

Reset Selection to Defaults

- Tower Selection
- Solvent Prewash 1
- Solvent Prewash 2
- Dispense Wash
- Dispense Pumps
- Dispense Settings
- Solvent Postwash 1
- Solvent Postwash 2

**Dispense Settings**

Draw Speed (µL/min): 1250

Dispense Speed (µL/min): 3000

Draw Needle Depth Offset (mm): 0.0

Use Needle Depth Offset for Dispense:

Viscosity Delay (s): 3

Air Gap (% Syringe Size): 0

Overfill (% Syringe Size): 5

Slower draw speeds reduce cavitation bubbles

Check box to enable washes or pumps, then select repetitions.  
"Tower Selection" and "Dispense Settings" are always enabled.

Reset Selection to Defaults

- Tower Selection
- Solvent Prewash 1
- Solvent Prewash 2
- Dispense Wash
- Dispense Pumps
- Dispense Settings
- Solvent Postwash 1
- Solvent Postwash 2

**Dispense Settings**

Draw Speed (µL/min): 1250

Dispense Speed (µL/min): 3000

Draw Needle Depth Offset (mm): 0.0

Use Needle Depth Offset for Dispense:

Viscosity Delay (s): 3

Air Gap (% Syringe Size): 0


Overfill (% Syringe Size): 5

Overfill keeps bubbles in the syringe and out of the vial



# Eliminating Syringe Bubbles Caused by In-Vial Vacuum

- Pre-slit septa in vial caps
- New needle guide design
- Extended nose on needle guide opens pre-slit septa
- Relieves vacuum inside vial
- Split reseals after needle withdraw
- Needle guide does not affect standard septa
- Not recommended for highly volatile samples



# Elimination of Syringe Bubble

- Pre-slit septa and new needle guide
- Slower draw speed (1250 uL/min)
- 5% Syringe overfill



# Running WorkBench Methods Using the Sequence Queue

The screenshot displays the WorkBench software interface for 'Instrument 1 (online): Method and Run Control GCC\_4.M, GCC\_2.s'. The main window is titled 'Method and Run Control' and shows the instrument is 'Ready'. The 'Instrument Status' panel on the left shows a 24°C temperature and 'Off' status. The 'Sequence Queue' panel on the right shows the active queue with two sequences: 'GCC\_2' (Pending) and 'GCC\_4.M' (Pending). The 'History Queue' panel shows a list of executed sequences with their completion times and statuses. The 'Current Logbook File INSTR1.LOG' panel at the bottom shows a table of event messages.

Name	Time entered into Qu...	Estimated Completion...	Status
GCC_2	9/21/2011 7:33:20 PM		Pending

Name	Time completed	Status
GCC_2	9/21/2011 7:33:04 PM	Completed
GCC_2	9/21/2011 7:30:27 PM	Completed
GCC_2	9/21/2011 7:24:16 PM	Completed
GCC_2	9/21/2011 7:17:46 PM	Aborted
GCC_2	9/21/2011 7:14:52 PM	Aborted
GCC_2	9/21/2011 7:13:28 PM	Aborted
GCC_2	9/21/2011 7:07:52 PM	Aborted

Module	#	Event Message	Date	Time
CP Command		Queued item GCC_2 changed state to Completed.	09/21/11	19:33:04
Method		Loading Method GCC_4.M	09/21/11	19:32:57
Sequence		GCC_2.s completed	09/21/11	19:32:54
Logbook		Sequence Logbook file saved in GCC_2.LOG	09/21/11	19:32:54
Sequence		GCC_2.s started	09/21/11	19:31:09
CP Command		Queued item GCC_2 changed state to Running.	09/21/11	19:31:09
CP Macro		Saving Sequence GCC_2.s	09/21/11	19:31:08
CP Command		Queued item GCC_2 changed state to Submitted.	09/21/11	19:31:07

- Methods are run using Sequence Queue
- Select the Method
- Choose the number of sample vials to run
- WorkBench will keep track of resource usage during sequence

# Method IP585 – 10 mg/kg Standard Preparation

Steps	WorkBench Prep of 10 mg/kg IP585 Standard
1	Solvent wash 250 uL syringe (front tower)
2	Add 90 uL dodecane to vial
3	Solvent wash 25 uL syringe
4	Add 10 uL working standard solution (100 mg/kg)
5	Solvent wash 25 uL syringe (back tower)
6	Add 1 uL ISTD solution (C17:0d33)
7	Mix



# Method IP585 – 10 mg/kg Standard Preparation

Steps	WorkBench Prep of 10 mg/kg IP585 Standard
1	Solvent wash 250 uL syringe (front tower)
2	Add 90 uL dodecane to vial
3	Solvent wash 25 uL syringe
4	Add 10 uL working standard solution (100 mg/kg)
5	Solvent wash 25 uL syringe (back tower)
6	Add 1 uL ISTD solution (C17:0d33)
7	Mix



# Method IP585 – 10 mg/kg Standard Preparation

Steps	WorkBench Prep of 10 mg/kg IP585 Standard
1	Solvent wash 250 uL syringe (front tower)
2	Add 90 uL dodecane to vial
3	Solvent wash 25 uL syringe
4	Add 10 uL working standard solution (100 mg/kg)
5	Solvent wash 25 uL syringe (back tower)
6	Add 1 uL ISTD solution (C17:0d33)
7	Mix




# Method IP585 – Complete Automation of Calibration Standards and Sample Preparation


- 2 Sets of Five Calibration Standards
  - Low Level Set: 0, 2, 46, 8, and 10 mg/kg FAMEs in C12
  - High Level Set: 20, 40, 60 , 80 and 100 mg/kg FAMEs in C12
- 3 Jet Fuel Samples with Different Amounts of Total FAME
  - 3-Step WorkBench Sample Prep:
    - Add 100 uL sample to empty vial
    - Add 1 uL ISTD to sample
    - Mix
  - Each sample prepped in duplicate to measure repeatability (r)
  - One injection of each sample on 5975C GC/MS system
    - Configured for IP method 585

# Comparison of Manual and WorkBench Calibration

Methyl Stearate – Manual Prep  
using 1 mL Cal Stds




Methyl Stearate – 7696A Prep  
Using 100 uL Cal Stds



**Better calibrations using fewer amounts of expensive compounds**



# WorkBench Prep Sample Precision



# WorkBench Prep Sample Precision

Consistent results obtained for each sample

Greatly exceeds IP585 repeatability (r) specification at all three levels of FAME content

Jet Fuel Sample A

	Run 1	Run 2	Avg	Std Dev	RSD	r	r (IP585)
C16:0	17.53	17.25					
C17:0	16.71	16.39					
C18:0	16.56	16.32					
C18:1	16.93	16.70					
C18:2	17.22	16.99					
C18:3	17.72	17.50					
<b>Total</b>	<b>102.67</b>	<b>101.15</b>	<b>101.91</b>	<b>1.07</b>	<b>1.05%</b>	<b>1.52</b>	<b>19.63</b>

Jet Fuel Sample B

	Run 1	Run 2	Avg	Std Dev	RSD	r	r (IP585)
C16:0	6.05	6.17					
C17:0	0.00	0.00					
C18:0	2.10	2.13					
C18:1	9.69	9.85					
C18:2	28.52	28.90					
C18:3	4.67	4.71					
<b>Total</b>	<b>51.03</b>	<b>51.76</b>	<b>51.40</b>	<b>0.52</b>	<b>1.00%</b>	<b>0.73</b>	<b>11.39</b>

Jet Fuel Sample C

	Run 1	Run 2	Avg	Std Dev	RSD	r	r (IP585)
C16:0	1.20	1.30					
C17:0	0.00	0.00					
C18:0	0.20	0.20					
C18:1	0.90	1.10					
C18:2	1.90	2.00					
C18:3	0.90	0.90					
<b>Total</b>	<b>5.10</b>	<b>5.50</b>	<b>5.30</b>	<b>0.28</b>	<b>5.34%</b>	<b>0.40</b>	<b>3.86</b>



# Summary

- Six simple steps used to translate a manual sample prep to an automated WorkBench Prep
- WorkBench hardware provide tools need for many liquid sample prep techniques
  - 2 syringe pipettes, wash solvents, waste containers, heating and cooling stations, mixing station
- WorkBench software provides easy, graphical user interface to build and run automated sample prep techniques
  - resource layout, drag-and-drop functions, Sequence Queue, logbooks
  - software automatically keeps track of resource usage
    - Warns user when resource are low
    - Warns user when resources are full (i.e. waste vials)