



7693A Automatic Liquid Sampler Supplies Quick Reference Guide

Agilent Supplies for Agilent Instruments

Agilent Technologies is committed to optimizing your laboratory's productivity. As part of this long-term commitment, we have designed our supplies to match our instruments precisely. You get consistent results and optimum instrument performance, all day, every day.

For your convenience, the reference list below contains the most commonly ordered supplies and parts.

Description	Part No.	Notes
7693A Replacement Parts and Supplies		
Gripper finger caps, 16/pk	G4514-60710	Need occasional replacement; keep a spare set on hand.
Injector mounting post	G4513-20561	Needed if one injector is to be moved between the front inlet and the rear inlet.
Dual parking post for autosampler	05890-61525	Provides a convenient place to park your injector while performing inlet maintenance.
Needle support insert, standard	G4513-40525	When working frequently with dirty samples, replace this part regularly to maintain smooth syringe movement.
Needle support insert, on-column	G4513-40529	
7693A vial trays with colored labels, 3/pk with 12 label tags (4 colors)	G4525A	Extra set is very useful for working on the bench and then exchanging with racks on autosampler
Vial rack label kit	G4525-60701	
Vial rack labels kit, red, 3/pk	G4525-60702	
Vial rack labels kit, yellow, 3/pk	G4525-60703	
Vial rack labels kit, green, 3/pk	G4525-60704	
ALS main cable assembly	G4514-60610	Provides added convenience in moving injector between front and rear inlets; eliminates unplugging and replugging cable in back of instrument.
Agilent Blue Line Autosampler Syringes for 7693A ALS		
Syringe, 1 µL, fixed needle, 23/42/Cone	G4513-80229	Plunger in needle design provides accuracy for small volume injections.
Syringe, 5 µL, fixed needle, 23-26s/42/HP	G4513-80206	Fitted plunger, fixed tapered needle.
Syringe, 5 µL, on-column, removable needle, barrel only	G4513-80207	For on-column injections; requires needle.
Replacement needle, 0.32 mm on-column, 5 µL syringe, 3/pk	G4513-80214	
Replacement needle, 0.25 mm on-column, 5 µL syringe, 3/pk	G4513-80217	

Continued on the next page.



Agilent Technologies

Description	Part No.	Notes
Agilent Blue Line Autosampler Syringes for 7693A ALS		
Syringe, 10 µL, PTFE tip, fixed needle, 23-26/42/HP	G4513-80203	Recommended for best precision and accuracy; PTFE-tipped plunger gives lowest carryover.
Syringe, 10 µL, fixed needle, 23-26s/42/HP	G4513-80204	For general use; fitted plunger and tapered needle
Syringe, 10 µL, PTFE tip, removable needle, 23-26/42/HP	G4513-80233	Both plunger and needle can be replaced for extended use.
Replacement needle for 10 µL syringe, 23-26/42/HP, 3/pk	G4513-80231	
Syringe, 25 µL, PTFE tip, fixed needle, 23/42/HP	G4513-80228	
Syringe, 50 µL, PTFE tip, fixed needle, 23/42/HP	G4513-80221	
Syringe, 100 µL PTFE tip, fixed needle, 23-26s/42/HP	G4513-80222	
Agilent Blue Line Syringes for Enhanced Sample Handling Carriage (G4521A)		
Syringe, 250 µL, PTFE tip, fixed needle, 23/42/HP	G4513-60560	250 µL size is more precise at lower volumes and more versatile.
Syringe, 500 µL, PTFE tip, fixed needle, 23/42/HP	G4513-60561	500 µL size is more popular; dispenses solvent quicker.
For a complete listing of syringes, replacement needles, and replacement PTFE tip plungers, go to www.agilent.com/chem/syringes		
2 mL Wide Opening Snap Top Glass Vials		
Pre-assembled clear screw top vial convenience pack, 100/pk Includes clear screw top vials, blue screw caps, PTFE/red silicone septa	5182-0553	
Pre-assembled amber screw top vial convenience pack, 100/pk Includes amber screw top vials with write-on spot, green screw caps, PTFE/red silicone septa	5182-0554	Amber for light-sensitive samples.
Pre-assembled clear screw top vial convenience pack, 100/pk Includes clear screw top vials, blue screw caps, PTFE/silicone septa	5182-0557	
Non-assembled screw top vial convenience pack, 500/pk Includes clear screw top vials, blue screw caps, PTFE/silicone septa in a 6 drawer storage box	5182-0734	
Non-assembled screw top vial convenience pack, 500/pk Includes clear screw top vials with write-on spot, blue screw caps, PTFE/red silicone rubber septa in 6 drawer storage box	5182-0867	
Clear screw top vials, 100/pk	5182-0714	
Amber screw top vials with write-on spot, 100/pk	5182-0716	Amber for light-sensitive samples.
Blue screw caps, PTFE/red silicone septa, 100/pk	5182-0717	
Clear, silanized screw top vials, 100/pk	5183-2070	Silanized for active or "sticky" compounds.
2 mL Wide Opening Crimp Top Glass Vials		
Non-assembled crimp top convenience pack, 500/pk Includes clear crimp top vials, silver aluminum crimp caps, PTFE/red rubber septa in 6 drawer storage box	5181-3400	
Clear crimp top vials, 100/pk	5181-3375	
Amber crimp top vials with write-on spot, 100/pk	5181-3376	Amber for light-sensitive samples.
Silver-tone aluminum crimp caps, PTFE/red rubber septa, 100/pk	5181-1210	
Crimping and Decapping Tools		
Electronic crimper, 11 mm, with 4.8v rechargeable battery pack and charger	5062-0207	Ideal for preparing large runs of samples in crimp vials.
Electronic decapper, 11 mm, with 4.8v rechargeable battery pack and charger	5062-0209	Ideal for sample recovery or recycling of large quantities of vials.
Manual crimper, 11 mm	8710-0979	
Manual decapper, 11 mm	5181-1213	

Continued on the next page.

Description	Part No.	Notes
2 mL Wide Opening Snap Top Glass Vials		
Non-assembled snap top vial convenience pack, 500/pk Includes clear snap top vials, clear polypropylene snap caps, PTFE/red silicone septa in a 6 drawer storage box	5182-0547	Crimp caps may also be used on snap cap vials.
Clear crimp/snap top vials, 100/pk	5182-0544	
Amber crimp/snap top vials with write-on spot, 100/pk	5182-0545	Amber for light-sensitive samples.
Green polyethylene caps, clear PTFE/red silicone septa, 100/pk	5182-3457	
Microvolume Inserts for Wide Opening Vials		
250 µL glass inserts with polymer feet and mandrel interior, 100/pk	5181-1270	
250 µL pulled-point glass inserts, 100/pk	5183-2085	
400 µL glass flat bottom inserts, 500/pk	5181-3377	
Shell Vials		
Clear, glass shell vials and caps, 200/pk	5190-1570	New economical choice for use with 7693A autosamplers.
Wash Vials (also for standards, diluents)		
4 mL wash vials with fill markings and caps, 25/pk	5182-0551	Autosampler holds 15 wash vials, but unit is shipped with a starter pack of 10.
Diffusion caps for 4 mL vials, 12/pk	07673-40180	Eliminates septum coring with repeated entry.
Septa for 4 mL vials, 144/pk	9301-1031	
Vial Accessories		
Rack for 12 mm, 2 mL vials, 5/pk	9301-0722	Holds 50 vials
Barcode labels numbered 1 to 1000, 1000/pk	5958-9450	
Agilent offers an extensive selection of vials, caps and septa in various types, materials and colors. For a complete listing, go to www.agilent.com/chem/vials		

**To learn more about Agilent's complete offering
of parts and supplies for your 7693A ALS:**

Visit www.agilent.com/chem/7693Asupplies

Contact your local Agilent Representative

Contact your local Agilent Authorized Distributor

Visit www.agilent.com/chem/contactus
to find an Agilent customer center in your country

Agilent Technologies shall not be liable for errors contained
herein or for incidental or consequential damages in
connection with the furnishing, performance, or use of this
material. Information, description, and specifications in this
publication are subject to change without notice.

© Agilent Technologies, Inc. 2009
Printed in U.S.A. August, 2009
5990-4586EN



Agilent Technologies