

AGILENT ASSAYMAP SOLUTION ASSAYMAP SYRINGE REPLACEMENT KIT

Maintain Optimal AssayMAP
Bravo Performance

Highlight

- Rapid syringe replacement to minimize down time.



Rapid AssayMAP syringe replacement

The AssayMAP Syringe Replacement kit includes the tools required to replace an AssayMAP syringe, step-by-step instructions, and 10 replacement syringes. The procedure can be performed rapidly by someone who is familiar with the AssayMAP Bravo Platform and the Protein Sample Prep Workbench. After replacing a syringe, perform the Agilent AssayMAP Syringe Test Kit (G5496-60050) to verify the optimal performance of the AssayMAP Bravo.

Product description	Part number	Kit components
Syringe Replacement Kit	G5409-68002	All tools required to replace a syringe 10 pack of AssayMAP syringes Detailed step-by-step instructions

Required by not included in the kit	Part number	Description
Reservoir Plates	Seahorse 201280-100 Seahorse 201282-100 Agilent G5498B#049	Labware options are listed. A single reservoir plate is required.

www.agilent.com/lifesciences/AssayMAP

This item is intended for Research Use Only.
Not for use in diagnostic procedures. This
information is subject to change without notice.

© Agilent Technologies, Inc. 2015
Published in the USA, August 14, 2015
5991-6148EN

