


Make Your Mass Spectrometry More Powerful

Special upgrade promotion:
Agilent MassHunter software



There's never been a better time to upgrade to the latest MassHunter software for GC/MS and LC/MS systems

Here's what's new:


- Agilent MassHunter Acquisition for LC/TOF and Q-TOF, Quantitative Analysis, and BioConfirm now contain technical controls to support 21 CFR Part 11 and Annex 11 compliance.
- A common Control Panel is shared across the three programs that allows the setup of user roles, user permissions, audit trails, checksums, and secure content management across all software modules. These advanced security features and built-in technical controls minimize security risks while helping to preserve the integrity of the data.
- MassHunter Acquisition for LC/TOF and Q-TOF now fully controls Agilent's 2D-LC system.
- MassHunter Acquisition for 6400 Series LC/TQ includes compliance feature enhancements including support for ECM 3.5/3.6.
- MassHunter for GC/MS systems now supports the 8860 and the innovative 8890 GC.
- GC/MS MassHunter includes compliance feature enhancements including support for ECM 3.5/3.6, support for QuickProbe for GC/MSD systems, and CI source support for GC/Q-TOF.

Discounts available for Agilent MassHunter software customers

This is a limited-time offer, so claim your discount now.

- Save 45% to 60% on MassHunter software products.
- Save 50% on MassHunter-ready computer bundles (Windows 10 Professional included).

To claim your discount visit
www.agilent.com/chem/masshunter_upgrade



Some of MassHunter products eligible for discount:*

- G1668AA MassHunter-ready PC workstation for MassHunter-based LC/MS systems
- G6842AA MassHunter GC/MS software upgrade bundle
- G1701FA GCMSD MassHunter SW Upgrade
- M5961AA Upgrade MassHunter Workstation for LC/TOF and LC/Q-TOF

* For a complete list of eligible products, visit www.agilent.com/chem/masshunter_upgrade

Additional reasons to upgrade:

- Windows 7 support ended in January 2020; you can upgrade your operating system with MassHunter-ready computer bundles, complete with Windows 10 Professional.
- You'll receive a one year software maintenance agreement (SMA) with the purchase of the upgrade products
- You'll receive free MRM database updates for three years from date of purchase, allowing you to keep your experiments current with access to the latest content and updated methods.

To claim your discount visit www.agilent.com/chem/masshunter_upgrade

DE44316.8855671296

This information is subject to change without notice.

© Agilent Technologies, Inc. 2019, 2021
Published in the USA, May 6, 2021
5994-0914EN

