



Save Time and Money with the Agilent JetClean Self-Cleaning Ion Source for GC/MS



Save up to 90%*
on instrument downtime

JetClean self-cleaning system
maximizes sample throughput



Save time
by avoiding reruns

A clean ion source ensures
data consistency




Save >75%
on cleaning costs

JetClean automated cleaning
requires virtually no user intervention

* Saving rate depends on the application.

Free yourself from manual ion source cleaning and lab downtime

The patented Agilent JetClean self-cleaning ion source greatly reduces—or even eliminates—the need for source cleaning on Agilent single and triple quadrupole GC/MS systems. A carefully controlled flow of hydrogen ensures that the ion source is kept free from contamination, protecting the integrity of your results.



Calculate how much your lab could save

Download our e-book and use the built-in cost calculator

For more information or a quote, visit www.agilent.com/chem/jetclean

This information is subject to change without notice.