



Agilent 7010B and 7000D Triple Quadrupole GC/MS Supplies

Agilent supplies for Agilent instruments

Agilent Technologies is committed to optimizing your laboratory's productivity, so we have produced this list of the most commonly ordered supplies and parts for the 7010B and 7000D Triple Quadrupole GC/MS. Keep this list handy so you can quickly find the supplies you need and minimize instrument downtime.


Foreline pump supplies	Part No.
Foreline pump oil, P3, 1 qt	6040-0834
Rough pump inlet flange	0905-1463
Oil drain kit for mist filter	3162-1057
Pump oil drip pan	G1946-00034
Oil mist filter for RV5 pump	G6600-80043
Quiet cover for RV series pumps	G3199C

Gas Clean filters	Part No.
Gas Clean purifier, carrier gas	CP17973
Gas Clean purifier kit (includes one connecting unit 1/8 in and two GC/MS filters)	CP17976
Gas Clean purifier kit (includes one connecting unit 1/4 in and two GC/MS filters)	CP17977
Gas Clean purifier installation kit	CP17978

Gas filters	Part No.
Chemical ionization gas purifier	5190-9071
Big universal trap, 1/8 in fittings, helium (Ar/Me)	RMSH-2
Big universal trap, 1/8 in fittings, nitrogen	RMSN-2
Big universal trap, 1/8 in fitting, hydrogen	RMSHY-2
Mounting clip	UMC-2

Maintenance supplies	Part No.
Abrasive sheets	5061-5896
Abrasive powder, aluminum oxide	393706201
Cloths, lint-free, 23 x 23 cm, 100% cotton, 15/pk	05980-60051
Cloths, lint-free, 23 x 23 cm, 100% cotton, 300/pk	9310-4828
Cotton swabs	5080-5400
Nylon gloves, lint-free, large	8650-0030
Nylon gloves, lint-free, small	8650-0029
High vacuum grease, 25 g	6040-0289
Electron multiplier	G7002-80103
4 Turn EI filament	G7005-60061
CI filament – 2/pk	G7005-60072
Filament assembly for 7010B	G7002-60001
Micro-Ion Gauge E with flange	G1960-80303
Replacement glass bulb for PFTBA and PFDTD test sample	G3170-80002

Ferrules and connectors	Part No.
Self Tightening column nut	5190-6194
Self Tightening column nut, MS	5190-5233
Ferrule, 0.4 mm id, 15% graphite/85%Vespel, 0.1 to 0.25 mm column	5181-3323
Ferrule, 0.5 mm id, 15% graphite/85%Vespel, 0.32 mm column	5062-3514



Agilent supplies—to help you get the results you need

For greater confidence in your chromatography results, you can optimize your total system with quality Agilent columns and supplies. Our wide range of columns, supplies, kits, and accessories is designed, manufactured, and tested to rigorous Agilent specifications, under a quality system registered to ISO 9001. Why risk compromising your analytical results with anything less than genuine Agilent supplies?

When choosing your columns, Agilent J&W GC columns offer the lowest bleed levels, the best inertness for acids/bases/mixed functional compounds, and the tightest column-to-column reproducibility.

For Agilent's complete listing of columns and supplies, visit www.agilent.com/chem/supplies

This information is subject to change without notice.

© Agilent Technologies, Inc. 2017
Printed in the USA, December 7, 2017
5991-8719EN

 **Agilent**
Trusted Answers