

Agilent 5100/5110 시리즈 ICP-OES

Dual View ICP-OES 소모품

애질런트 테크놀로지스는 귀하의 실험실 생산성을 최적화하기 위해 최선을 다하고 있으며, 이에 Agilent 5100/5110 SVDV/VDV 및 RV ICP-OES 시스템에 가장 많이 쓰이는 소모품과 부품 리스트를 발행했습니다. 이 리스트를 보유하면 필요한 소모품을 빠르게 찾아볼 수 있고 기기 가동 중단 시간을 최소화할 수 있습니다.



Nebulizer 및 소모품	부품 번호
Concentric glass nebulizer, SeaSpray U Series with Ezylok, standard for 5100/5110.	G8010-60255
Concentric glass nebulizer, Conikal U Series with Ezylok Recommended for organics.	G8010-60270
OneNeb Series 2 inert concentric nebulizer. Suitable for most sample types including organic solvents, HF acid digests and high TDS samples.	G8010-60293
Sample inlet fitting for OneNeb Series 2 nebulizer. Includes capillary connector with captive sample capillary tubing 0.5 mm id, 1/16 inch od, 750 mm long.	G8010-60345
Cleaning tool for concentric nebulizer.	G3266-80020
Sample inlet (Unifit) fitting for concentric glass nebulizer 0.75 mm id x 700 mm L, 10/pk.	G8010-80035
Sample inlet (Unifit) fitting for concentric glass nebulizer 0.5 mm id x 700 mm L, 10/pk.	G3266-80012
Gas inlet connector (Ezylok) for concentric nebulizer, 1/pk.	G3266-80015
PFA tubing, nebulizer gas supply, per m.	0890-2603

펌프 및 연결 튜빙	부품 번호
Peristaltic pump tubing, PVC, white/white, 12/pk.	3710034400
Peristaltic pump tubing, PVC, blue/blue, 12/pk.	3710034600
Peristaltic pump tubing, PVC, black/black, 12/pk.	3710027200
Peristaltic pump tubing, PVC, orange/white, 12/pk.	3710046900
Y piece connector for online addition of internal standard or ionization buffer.	1610132400
Barb connector, 1/16" to 1/8", connects peristaltic pump tube for waste to drain tubing.	1610088600
Tubing, waste line from peri pump tube, per m.	3710024600
Tubing, waste line from spill tray, per m.	3710033400
Nebulizer capillary tubing, per m.	2410020500

스프레이 챔버	부품 번호
Spray chamber, double pass glass cyclonic, standard for 5100/5110.	G8010-60256
Spray chamber, single pass glass cyclonic.	G8010-60271
Inert double pass spray chamber, white, rec. for HF samples.	G8010-80034
Inert single pass spray chamber, white, rec. for HF samples.	G8000-70008
Drain connector for spray chamber waste (Unifit) 0.86 mm id x 2.0 mm od x 480 mm L, 3/pk.	G8010-80036
Seal for spray chamber nebulizer fitting (Helix), 4/pk.	G8010-80042
Locking screw for s/chamber neb. fitting (Helix).	G8010-80043

Advanced Valve System(AVS) 소모품	부품 번호
기타 AVS 소모품	
Rinse (20 L) and waste (10 L) containers for AVS 6/7 (1 each).	G8494-60005
Internal standard container (5 L) kit for AVS 7, 1/pk.	G8495-60005
Reagent uptake capillary tube 1.6 mm od x 1.0 mm id x 1400 mm long with heavy PTFE sinker, 1/pk.	G8494-60007
Check valve outlet for AVS switching valve, non-metallic, 1/pk.	G8010-60300
Mounting clamp for AVS 4/6/7 switching valve, 1/pk.	G8494-60004
AVS용 튜빙 키트	
Sample inlet connection kit for AVS switching valve, 1/pk.	G8010-60303
Nebulizer connection kit for AVS switching valve, 1/pk.	G8010-60304
Autosampler connection kit for AVS switching valve, 1/pk.	G8010-60306
Internal standard connection kit for AVS 7 switching valve, 1/pk.	G8010-60307
Pump outlet connection kit for AVS switching valve, 1/pk.	G8010-60308
Nebulizer adaptor kit for AVS switching valve, 1/pk.	G8010-60333
Waste tube kit for 4 port AVS switching valve, 1/pk.	G8010-60334
Sample/rinse connection kit for AVS switching valve, 2/pk.	G8010-60335
Autosampler peri-pump connection kit, 1/pk. Allows bypass of the switching valve when required e.g. for wavelength calibration.	G8010-60348
AVS 시료 루프	
AVS sample loop, 1.0 mL volume, 1/pk. (Standard loop supplied with the AVS).	G8010-60305
AVS sample loop, 0.25 mL volume, 1/pk.	G8010-60338
AVS sample loop, 0.5 mL volume, 1/pk.	G8010-60339
AVS sample loop, 1.5 mL volume, 1/pk.	G8010-60340
AVS sample loop, 2.0 mL volume, 1/pk.	G8010-60341
AVS sample loop, 2.5 mL volume, 1/pk.	G8010-60342
AVS sample loop, 3.0 mL volume, 1/pk.	G8010-60343
AVS sample loop, 3.5 mL volume, 1/pk.	G8010-60344
AVS 7/6/4 밸브 소모품	
Spare valve, 7 port, 2 position for AVS 7 switching valve, 1/pk.	G8495-60000
Stator and rotor for 7 port AVS 7 switching valve, 1/pk.	G8495-60003
Rotor seal for 6-7 port valve used with AVS 6/7 switching valve, 1/pk.	G8494-60002
Spare valve, 6 port, 2 position for AVS 6 switching valve, 1/pk.	G8494-60000
Stator and rotor for 6 port AVS 6 switching valve, 1/pk.	G8494-60003
Spare valve, 4 port, 2 position for AVS 4 switching valve, 1/pk.	G8493-60000
Stator and rotor for 4 port AVS 4 switching valve, 1/pk.	G8493-60003
Rotor seal for 4 port AVS 4 switching valve, 1/pk.	G8493-60002

토치 및 소모품	부품 번호
Axial 및 radial 뷰에 사용하는 5100/5110 SVDV/VDV 토치	
Easy-fit torch, one piece with 1.8 mm id injector and an extended outer tube with slot, standard for SVDV/VDV models.	G8010-60228
Easy-fit torch, demountable with 1.8 mm id injector, an extended outer tube with slot & replaceable inner/outer tube set.	G8010-60236
Easy-fit torch, demountable, high solids with 2.4 mm id injector, an ext. outer tube with slot & replaceable inner/outer tube set.	G8010-60232
Easy-fit torch, demountable, semi-volatile organics, with 1.4 mm id inj., ext. outer tube with slot & replaceable inner/outer tube set.	G8010-60233
Easy-fit torch, inert, demountable for HF digests, with 1.8 mm id ceramic injector, an ext. outer tube with slot and replaceable inner/outer tube set.	G8010-60231
Quartz outer tube set for demountable torch, comprises ext. length outer tube with slot, intermediate tube (integrated) and top seal.	G8010-60263
전용 radial 뷰를 위한 5100/5110 RV 토치	
Easy-fit torch, one piece with 1.4 mm id injector. Has a shorter length outer tube (no slot), standard for RV models.	G8012-60000
Easy-fit torch, demountable with 1.4 mm id injector, shorter outer tube (no slot) and replaceable inner/outer tube set.	G8010-60237
Easy-fit torch, demountable, volatile organics, with 0.8 mm id injector, shorter outer tube (no slot) and replaceable inner/outer tube set.	G8010-60235
Easy-fit torch, inert, demountable for HF digests, with 1.8 mm id ceramic injector, shorter outer tube (no slot) and replaceable inner/outer tube set.	G8010-60234
Quartz outer tube set for demountable torch, comprises shorter length outer tube (no slot), intermediate tube (integrated) and top seal.	G8010-60264
기타 5100/5110 토치 소모품	
Torch clamp, to clamp the spray chamber to the ball joint on the torch (not suitable for inert torch).	G8000-64140
Torch clamp for inert, demountable torch. Clamps the inert spray chamber to the base of the inert torch (not suitable for standard torch).	G8010-60227
Torch cleaning stand for 5100/5110 Series torches. Designed to make the torch cleaning process easier and safer, also preventing damage to the torch.	G8010-68021
Storage rack for 5100/5110 Series torches. Holds up to 5 torches, multiple racks can be stacked on top of each other if required.	G8010-67000

작동 소모품 키트	부품 번호
Includes supplies required to support up to 12 months routine operation. Extra supplies may be required for labs. that have high workload or high matrix samples.	
Operating supplies kit for 5100/5110 DV ICP-OES.	G8010-67001
Operating supplies kit for 5100/5110 RV ICP-OES.	G8012-67001
Each kit includes:	
<ul style="list-style-type: none"> • Demountable torch • 2 quartz outer tube sets • 4 pks sample + waste peristaltic pump tubing • Internal standard pump tubing and connectors • SeaSpray nebulizer (U series) • Unifit connectors (0.75 mm id) for nebulizer • EzyLok gas connector • Double pass glass cyclonic spray chamber • Unifit connectors for spray chamber drain • Spare clamp for the ball joint • Replacement pre-optic window (DV kit includes supplies for both radial and axial windows) • Barb connectors for waste line • Nebulizer capillary tubing • PFA tubing for nebulizer gas supply • Drain tubing for waste liquid from spray chamber • Drain tubing for sample tray and seal for spray chamber nebulizer fitting (Helix) 	

응용 키트	부품 번호
HF 분해용 비활성 시료 주입 키트	
Kit includes:	
<ul style="list-style-type: none"> • Inert double pass spray chamber with ball joint socket • OneNeb Series 2 nebulizer • Demountable, inert Easy-fit torch with ceramic injector (1.8 mm id) • Torch clamp for inert torch • SolvaFlex peristaltic pump tubes for sample (black/black tabs) and waste (gray/gray tabs) 	
Inert sample introduction kit for SVDV/VDV.	G8010-68007
Inert sample introduction kit for RV.	G8010-68006
Ketone 기반 유기물 시료 주입 키트	
Kit includes:	
<ul style="list-style-type: none"> • Double pass glass cyclonic spray chamber with ball joint socket • Conikal nebulizer • Demountable Easy-fit torch with 1.4 mm id injector • Solvent resistant drain tubing • Marprene peristaltic pump tubes for sample (white/white tabs) and waste (blue/blue tabs) 	
Ketone-based organics sample introduction kit for SVDV/VDV.	G8010-68009
Ketone-based organics sample introduction kit for RV.	G8010-68008
휘발성 유기물 시료 주입 키트	
Kit includes:	
<ul style="list-style-type: none"> • Conikal nebulizer • Demountable Easy-fit torch with 0.8 mm id injector • Solvent resistant drain tubing • Viton peristaltic pump tubes for sample (white/white tabs) and waste (blue/blue tabs) 	
Requires but does NOT include the IsoMist peltier cooled spray chamber.	
Volatile organics sample introduction kit for SVDV/VDV.	G8010-68005
Volatile organics sample introduction kit for RV.	G8010-68004
준휘발성 유기물 시료 주입 키트.	
Kit includes:	
<ul style="list-style-type: none"> • Double pass glass cyclonic spray chamber with ball joint socket • Conikal nebulizer • Demountable Easy-fit torch with 1.4 mm id injector • Solvent resistant drain tubing • PVC SolvaFlex peristaltic pump tubes for sample (black/black tabs) and waste (gray/gray tabs) 	
Semi-volatile organics sample introduction kit for SVDV/VDV.	G8010-68029
Semi-volatile organics sample introduction kit for RV.	G8010-68030

기타 소모품	부품 번호
Tubing clamp for sample tray. Secures the drain tubing to the sample tray.	G8010-40026
Argon gas inlet filter. Requires annual replacement. Customer replaceable.	G8010-60136
Replacement filter element for cooling air inlet. Suits the dusty environment air inlet filter option.	G8010-60189
Quartz window for the axial view pre-optic window.	G8010-68014
O-ring for the axial view pre-optic window.	G8010-68018
Pre-optics cone for the axial view pre-optic window.	G8010-68109
Esteel powder for cleaning pre-optics cone.	8110002800
Quartz window for the radial view pre-optic window.	G8010-68015
Argon humidifier accessory. Used to humidify the nebulizer gas for samples >3% TDS. Includes:	G8010-60346
• 2 x 0.6 m 4 mm od gas lines	
• Fitting plate and bracket	
• Installation instructions	

SPS 4 자동 시료 주입기 소모품	부품 번호
SPS 4 자동 시료 주입기 옵션	
High uptake volume rinse kit for SPS 4. Recommended for use when connected to switching valves. Allows rinse uptake rate to be increased to 25 mL/min.	G8410-68260
SPS 4 rinse tube kit, for low volatile organics. Allows use of low volatile organic solvents in rinse station.	G8410-68261
Pump drain kit for SPS 4 autosampler, 1/pk. Provides the components required to allow the liquid overflow to be pumped from the rinse station.	G8410-68262
예비 프로브	
SPS 4 probe 0.8 mm id, green marker, inert, FEP sleeved, 1/pk (Standard probe for ICP-OES).	G8410-80102
SPS 4 probe 1.0 mm id, blue marker, inert, FEP sleeved, 1/pk.	G8410-80103
SPS 4 probe 1.2 mm id, black marker, inert, FEP sleeved, 1/pk (required for viscous, high matrix and high TDS samples).	G8410-80104
행금 펌프용 연동 펌프 튜빙	
Peristaltic pump tubing, PVC SolvaFlex, 3 bridged, gray/gray tabs (12 tubes/pk). Standard tubing for the rinse pump, for use with aqueous or low volatile organic solvents.	3710049000
Peristaltic pump tubing, Santoprene, 3 bridged, gray/gray tabs (12 tubes/pk). For use with ketone-based solvents.	3710057500
Peristaltic pump tubing, Viton, 3 bridged, black/black tabs with flared ends (12 tubes/pk). For use with aromatic organic solvents.	3710069100
랙 및 튜브	
Standards rack 34 position (1/pk). Provides 12 places for 50 mL 30 mm od centrifuge tubes and 22 places for 18 mm od test tubes. (Standard configuration for ICP-OES).	G8410-40000
Sample rack for 13 mm od tubes, 90 positions (1/pk).	6610026400
Sample rack for 16 mm od tubes, 60 positions (1/pk).	6610025400
Sample rack for 20 mm od tubes, 40 positions (1/pk).	6610025500
Sample rack for 25 mm od tubes, 24 positions (1/pk).	6610026500
Sample rack for 30 mm od tubes (50 mL centrifuge tubes), 21 positions (1/pk).	6610026600
Rack overlay for 60-position sample rack (1/pk, recommend 2 per rack). Allows 13 mm and some 16 mm od tubes to be used in 60-position rack.	6610026100
Rack overlay for 40-position sample rack (1/pk, recommend 2 per rack). Allows 18 mm od tubes to be used in 40-position rack.	6610026000
Centrifuge tubes polypropylene, 50 mL (500/pk) screw-capped, 29 mm od.	190065200
Polypropylene tubes, 17 mm x 130 mm, 1000/pk.	00000001600L

VGA 77 소모품	부품 번호
Sodium tetrahydroborate reagent, 500 g.	8210029100
Reagent pump tubing, black/black, 12/pk.	3710027200
Sample pump tubing, purple/black, 12/pk.	3710027300
Tubing and connector kit.	9910061900
ICP-OES hydride module (recommended to dedicate a module to Hg determinations).	9910062200
Gas-liquid separator.	9910062000
Pump tube beds kit.	9910050700

MSIS 소모품	부품 번호
Reagent pump tubing, PVC, black/black, 12/pk.	3710027200
Waste pump tubing, PVC, black/white 12/pk.	3710068900
MSIS nebulizer bushing (with O-rings).	G8479-20020
O-rings for MSIS nebulizer bushing, 1 pair. Includes 1 pair of replacement O-rings for the nebulizer bushing (1 internal and 1 external O-ring).	G8479-67000
MSIS drain tube assembly, 1/pk. Includes: <ul style="list-style-type: none"> 0.2 m length of 1/8 inch id drain tubing with bushing for waste outlet Barb adapter Single peristaltic pump tube, 2 bridged PVC with black/white tabs 	G8479-60100
MSIS reaction cone assembly, 1/pk. Includes: <ul style="list-style-type: none"> Reaction cone Bushing Threaded barb adapter Tubing connector for sample inlet 	G8479-60300
MSIS reagent tubing kit. Includes: <ul style="list-style-type: none"> O-rings for nebulizer bushing (1 pair) Tubing connectors for reagent/sample inlets Capillary tubing 0.01 inch id x 0.03 inch od to join to peristaltic pump tubing for reductant inlet 	G8479-68100

단일 원소 ICP 표준물질	부품 번호
Supplied in 100 mL bottles, concentration 10,000 ppm (1%) in 5% nitric acid (unless indicated otherwise), 1 mL of solution contains 10 mg of metal.	
Aluminum ICP Standard.	5190-8352
Antimony ICP Standard in 1% HNO ₃ , 1% tartaric acid.	5190-8354
Arsenic ICP Standard.	5190-8356
Barium ICP Standard.	5190-8358
Beryllium ICP Standard.	5190-8360
Bismuth ICP Standard.	5190-8362
Boron ICP Standard in 1% NH ₄ OH.	5190-8364
Cadmium ICP Standard.	5190-8366
Calcium ICP Standard.	5190-8368
Cerium ICP Standard.	5190-8370
Cesium ICP Standard.	5190-8372
Chromium ICP Standard.	5190-8374
Cobalt ICP Standard.	5190-8376
Copper ICP Standard.	5190-8378
Dysprosium ICP Standard.	5190-8380
Erbium ICP Standard.	5190-8382
Europium ICP Standard.	5190-8384
Gadolinium ICP Standard.	5190-8386
Gallium ICP Standard in 5% HNO ₃ , 0.5% HCl.	5190-8388
Germanium ICP Standard in 5% HNO ₃ , trace HF.	5190-8390
Gold ICP Standard in 20% HCl.	5190-8392
Hafnium ICP Standard in 5% HCl.	5190-8394
Holmium ICP Standard.	5190-8396
Indium ICP Standard.	5190-8398
Iridium ICP Standard in 20% HCl.	5190-8400
Iron ICP Standard.	5190-8402
Lanthanum ICP Standard.	5190-8404
Lead ICP Standard.	5190-8406
Lithium ICP Standard.	5190-8408
Lutetium ICP Standard.	5190-8410
Magnesium ICP Standard.	5190-8412
Manganese ICP Standard.	5190-8414
Mercury ICP Standard.	5190-8416
Molybdenum ICP Standard in 1% NH ₄ OH.	5190-8418

단일 원소 ICP 표준물질(계속)	부품 번호
Neodymium ICP Standard.	5190-8420
Nickel ICP Standard.	5190-8422
Niobium ICP Standard in 2% HF.	5190-8424
Palladium ICP Standard in 10% HNO ₃ .	5190-8426
Phosphorus ICP Standard.	5190-8428
Platinum ICP Standard in 20% HCl.	5190-8430
Potassium ICP Standard.	5190-8432
Praseodymium ICP Standard.	5190-8434
Rhenium ICP Standard.	5190-8436
Rhodium ICP Standard in 20% HCl.	5190-8438
Rubidium ICP Standard.	5190-8440
Ruthenium ICP Standard in 20% HCl.	5190-8442
Samarium ICP Standard.	5190-8444
Scandium ICP Standard.	5190-8446
Selenium ICP Standard.	5190-8448
Silicon ICP Standard in H ₂ O.	5190-8450
Silver ICP Standard.	5190-8452
Sodium ICP Standard.	5190-8454
Strontium ICP Standard.	5190-8207
Sulfur ICP Standard in H ₂ O.	5190-8209
Tantalum ICP Standard in 2% HF.	5190-8211
Tellurium ICP Standard in 30% HCl.	5190-8213
Terbium ICP Standard.	5190-8215
Thallium ICP Standard.	5190-8217
Thulium ICP Standard.	5190-8219
Tin ICP Standard in 20% HCl.	5190-8221
Titanium ICP Standard in H ₂ O.	5190-8224
Tungsten ICP Standard in 5% HNO ₃ , trace HF.	5190-8226
Vanadium ICP Standard.	5190-8228
Ytterbium ICP Standard.	5190-8230
Yttrium ICP Standard.	5190-8232
Zinc ICP Standard.	5190-8234
Zirconium ICP Standard.	5190-8236

다원소 검량 표준물질	부품 번호
Calibration mix 1, 125 mL. Contains 100 mg/L Sb, Mo, Sn and Tl in 2% HNO ₃ + 0.5% HF.	6610030500
Calibration mix 2, 125 mL. Contains 100 mg/L Ag, Al, As, Ba, Be, Cd, Co, Cr, Cu, Mn, Ni, Pb, Se, Ti, Th, U, V and Zn in 5% HNO ₃ .	6610030600
Calibration mix majors, 125 mL. Contains 500 mg/L Ca, Fe, K, Mg and Na in 5% HNO ₃ .	6610030700
IntelliQuant 검량 용액	
IntelliQuant calibration kit. Used to update the pre supplied IntelliQuant calibration to achieve a greater level of accuracy for semi quantitative analysis. Includes multi-element standards No. 1, 2 and 3 plus a 5% HNO ₃ blank. NOTE: Hg standard is not part of this kit - if required, order part number 5190-8575.	5190-9425
IntelliQuant multi-element standard No.1. (Ready to use) 100 mg/L of 48 elements in 5% HNO ₃ (500 mL).	5190-9422
IntelliQuant multi-element standard No. 2. (Ready to use) 100 mg/L of Sb, Ge, Hf, Mo, Nb, Si, Ag, Ta, Sn, Ti, W, Zr; in 5% HNO ₃ with trace HF (500 mL).	5190-9423
IntelliQuant multi-element standard No. 3. (Ready to use) 100 mg/L Au, Ir, Os, Pd, Pt, Rh, Ru, Te; in 10% HCl (500 mL).	5190-9424
IntelliQuant single element standard for Hg. (Ready to use) 100 mg/L Hg, in 5% HNO ₃ (100 mL).	5190-8575
Calibration blank solution for ICP-OES. (Ready to use) Pure ASTM Type 1A water in 5% HNO ₃ (500 mL).	5190-7001

파장 검량 용액	부품 번호
ICP-OES wavelength calibration concentrate, 500 mL Contains 50 mg/L Al, As, Ba, Cd, Co, Cr, Cu, Mn, Mo, Ni, Pb, Se, Sr, Zn and 500 mg/L K in 5% HNO ₃ . Dilute x10 before use.	6610030000
ICP-OES wavelength calibration solution, 500 mL Contains 5 mg/L Al, As, Ba, Cd, Co, Cr, Cu, Mn, Mo, Ni, Pb, Se, Sr, Zn and 50 mg/L K in 5% HNO ₃ . Ready to use.	6610030100
ICP-OES calibration blank solution, 500 mL Pure ASTM Type 1A water with 5% HNO ₃ For use as cal. blank or diluent for the wavelength cal. concentrate.	5190-7001

내부 표준물질/이온화 억제제	부품 번호
Cesium Nitrate Buffer Solution, 1%. 10,000 ppm Cs in 5% nitric acid, 100 mL bottle.	5190-8343
Cesium Buffer Solution, 10%. 100,000 ppm Cs in 5% nitric acid, 500 mL bottle.	5190-9421
Yttrium Solution, 1%. 10,000 ppm Y in 5% nitric acid, 500 mL bottle. Also used for instrument optimization.	5190-8233
ICP Internal Standard, 125 mL bottle. 100 mg/L ⁶ Li, Sc, Y, In, Tb, Bi in 5% nitric acid.	6610030400

환경 표준물질	부품 번호
INTF-A Quality Control Standard. Interference check standard for preparation of US EPA CLP standards ICESA and ICSAB. Contains Al, Ca, Mg at 5000 ppm, Fe at 2000 ppm in 20% HCl, 500 mL.	5190-8599
ICV Quality Control Standard in 5% HNO ₃ with trace tartaric acid, 500 mL. Initial/continuing calibration verification standard for US EPA CLP analyses. Contains Ca, Mg, K, Na at 500 ppm, Al, Ba at 200 ppm, Fe at 100 ppm, Sb at 60 ppm, Co, V at 50 ppm, Ni at 40 ppm, Cu at 25 ppm, Zn at 20 ppm, Mn at 15 ppm, As, Cr, Ag, Tl at 10 ppm, Be, Cd, Se at 5 ppm, and Pb at 3 ppm.	5190-9408
QCSTD-27 Quality Control Standard, 100 mL. Quality control standard for environmental analyses. Contains Al, Sb, As, Ba, Be, B, Cd, Ca, Cr, Co, Cu, Fe, Pb, Mg, Mn, Mo, Ni, K, Se, Si, Ag, Sr, Na, Tl, Ti, V, Zn at 100 ppm.	5190-9418

애질런트는 또한 원자 분광기용 다양한 농도와 부피의 단일 원소 표준물질, 광범위한 다원소 표준물질, 이온화 버퍼, 금속 및 바이오디젤 표준물질, 용매 및 환경 표준물질을 제공합니다. NIST 추적성을 보장한 동시에, 애질런트의 인증 표준물질 및 표준 혼합물은 정확하고 신뢰할 수 있는 높은 가치의 검량을 제공합니다.

애질런트 부품 및 소모품 - 귀하가 원하는 결과를 얻을 수 있도록 도와드립니다

보다 확실한 분광기 결과를 위해 고품질의 애질런트 소모품으로 귀하의 전체 시스템을 최적화할 수 있습니다. 당사의 다양한 소모품, 키트, 액세서리 등은 ISO 9001의 인증을 받은 품질 관리 시스템에서 애질런트의 엄격한 기준에 맞춰 설계, 생산, 테스트됩니다.

왜 분석 결과에 영향주는 위험을 감수하며 정품 애질런트 소모품보다 품질이 낮은 제품을 선택하십니까?

AA, ICP-OES, ICP-MS 및 MP-AES용 중공 음극관 램프와 정품 애질런트 부품 및 소모품을 포함한 애질런트의 전체 분광기 소모품 목록을 살펴보시려면 다음 사이트를 방문하시기 바랍니다.

www.agilent.com/chem/atomic

www.agilent.com/chem

이 정보는 사전 고지 없이 변경될 수 있습니다.

© Agilent Technologies, Inc. 2019
2019년 6월 1일, 한국에서 인쇄
5991-4950KO

서울시 용산구 한남대로 98, 일신빌딩 4층 우)04418
한국애질런트테크놀로지스(주) 생명과학/화학분석 사업부
고객지원센터 080-004-5090 www.agilent.co.kr

